

# INSTRUCTION MANUAL

# **DNA Clean & Concentrator™-500**

Catalog Nos. **D4031 & D4032** 

# **Highlights**

- Simple, rapid recovery of ultra-pure DNA from large-scale sample sources.
- Unique column construction allows sample washing to be performed using a centrifuge or vacuum source.
- Eluted DNA is ideal for PCR, DNA sequencing, DNA transfection, DNA ligation, endonuclease digestion, RNA transcription, radiolabeling, etc.

# **Contents**

Product Contents	1
Specifications	1
Product Description	2-3
Buffer Preparation	4
Protocol	4
Troubleshooting	5
Ordering Information	6
List of Related Products	7-8

Page 1	
For Research Use Only	Ver. 1.2.0

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

### **Product Contents**

DNA Clean & Concentrator™-500 (Kit Size)	<b>D4031</b> (10 Preps.)	<b>D4032</b> (20 Preps.)	Storage Temperature
DNA Binding Buffer	100 ml	2 x 100 ml	Room Temp.
DNA Wash Buffer <sup>1</sup>	24 ml	48 ml	Room Temp.
DNA Elution Buffer	2 x 16 ml	2 x 50 ml	Room Temp.
Zymo-Spin™ VI Columns	10	20	Room Temp.
Instruction Manual	1	1	-

Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

# **Specifications**

- **DNA Purity** High-quality DNA ( $A_{260}/A_{280}$  >1.8) ideal for ligation, sequencing, labeling, PCR, microarray, transfection, transformation, and restriction digestion procedures.
- DNA Size Limits From ~50 bp to 23 kb.
- DNA Recovery Typically, up to 500 µg total DNA per column can be eluted into as little as 2 ml of low salt DNA Elution Buffer or water. For DNA 50 bp to 10 kb, the recovery is 70-90%. For DNA 11 kb to 23 kb, the recovery is 50-70%.
- **Sample Sources** DNA from enzymatic reactions (e.g., PCR, restriction endonuclease digestions), plasmid preparations, and impure preparations.
- Product Detergent Tolerance ≤ 5% Triton X-100, ≤ 5% Tween-20, ≤ 5% Sarkosyl,
   ≤ 0.1% SDS.
- Equipment Needed Swinging bucket centrifuge or vacuum source (optional).

Note: ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

<sup>&</sup>lt;sup>1</sup> Ethanol must be added prior to use as indicated on the **DNA Wash Buffer** label.

## **Product Description**

The <u>DNA Clean & Concentrator™-500</u> (DCC™-500) is designed for the rapid, large format purification and concentration of up to 500 µg of high quality DNA from samples including large-scale restriction endonuclease digestions and impure DNA preparations. The DCC™-500 features a single-buffer system that allows for efficient DNA adsorption onto the matrix of the supplied Zymo-Spin™ VI column. Simply add the specially formulated DNA Binding Buffer to your samples and transfer the mixtures to the supplied Zymo-Spin™ VI Column. The DNA is washed then eluted with a small volume (≥ 2 ml) of water or supplied DNA Elution Buffer using a centrifuge. The purified DNA is ideal for DNA ligation, sequencing, labeling, PCR, microarray, transfection, transformation, and restriction digestion procedures. The entire DNA purification/concentration procedure typically takes less than 20 minutes.

# Loading and washing the Zymo-Spin™ VI Column can be performed using any combination of the following:



Centrifuge: Zymo-Spin™ VI Column inside a 50 ml conical tube.




Vacuum: Zymo-Spin VI™ Column connected to a vacuum manifold.

#### Elute DNA Using a Centrifuge



DCC-500™ procedure.



Ultra-pure DNA from the **DNA Clean & Concentrator™** can be easily ligated for efficient cloning. Lanes: A: 3.5 kb bluntended vector; B: 7.4 kb blunt-ended DNA insert; C: Blunt-end ligation of A+B mixture.

#### **Available Formats**

	DCC™-5	DCC™-25	DCC™-100	DCC™-500	Genomic DCC™	ZR-96 DCC™-5
Name	Zymo-Spin™ I & IC	Zymo-Spin™ II & IIC	Zymo-Spin™ V	Zymo-Spin™ VI	Zymo-Spin ™ IC-XL	Zymo-Spin™ I-96
Capacity	5 μg/ prep.	25 μg/ prep.	100 μg/ prep.	500 μg/ prep.	10 μg/ prep.	5 μg/ prep.
Elution Vol.	≥ 6 µl	≥ 25 µl	≥ 150 µl	≥ 2 ml	≥ 10 µl	≥ 10 µl
Cat. Nos.	D4003, D4013	D4005, D4033	D4029, D4030	D4031, D4032	D4010, D4011	D4023, D4024

#### Typical DCC™ Applications

Post-PCR DNA Clean-up	Efficient desalting of DNA with the removal of DNA polymerases, primers and free dNTPs.
DNA Clean-up From Enzymatic Reactions	Efficient desalting of DNA with the removal of modifying enzymes, RNA polymerases, ligases, kinases, nucleases, phosphatases, endonucleases, etc.
Post-Reverse Transcription (RT) & cDNA Clean-up	Efficiently purifies DNA following RT, either as a DNA/RNA complex or as single stranded cDNA following chemical hydrolysis of the RNA template.
Plasmid DNA Clean-up	Efficiently purifies plasmid DNA from "home-made" preparations of cell free lysates or from commercial kits. Plasmid DNA purified and concentrated using the <b>DCC™</b> has proven an excellent substrate for high quality DNA sequencing.
Isotope and Dye Removal	Efficiently removes unincorporated fluorescent ( <i>i.e.</i> , AMCA, FITC, BIO, DIG, Cy3, Cy5, FAM, <i>etc.</i> ) and radiolabeled dNTP derivatives from DNA following <i>in vitro</i> labeling reactions.
Purification of M13 ssDNA	The <b>DCC</b> <sup>™</sup> can be used for the rapid isolation of single stranded M13 phage DNA directly from phage-infected <i>E. coli</i> culture supernatant.

- ✓ For purification of short DNA or RNA oligonucleotides ≥16 nt, use the Oligo Clean & Concentrator (D4060, D4061).
- ✓ For ChIP (Chromatin Immunoprecipitation) sample cleanup, use the ChIP DNA Clean & Concentrator (D5201, D5205) for high quality DNA from any step in a standard ChIP protocol.
- ✓ For post-cycle sequencing samples, use the ZR Sequencing DNA Clean-up Kit (D4050, D4051) for dye blob elimination.
- ✓ For samples containing PCR inhibitors, use the OneStep™ PCR Inhibitor Removal Kit (D6030, D6035).

#### **Selected Citations**

Li, N (2010) Whole genome DNA methylation analysis based on high throughput sequencing technology. *Methods*, *52* (3), 221-232. Khoo, S (2011) Acquiring genome-wide gene expression profiles in Guthrie card blood spots using microarrays. *Pathology*, *61* (1), 1-6. Lee, EJ (2011) Targeted bisulfite sequencing by solution selection and massively parallel sequencing. *Nucleic Acids Research*, *39*(19), e127, doi:10.1093/nar/gkr598

Papageorgiou, EA. (2009) Sites of differential DNA methylation between placenta and peripheral blood. *Am J Pathol, 174* (5), 1609-1618. Ferguson, AA et al. (2009) Retrofitting ampicillin resistant vectors by recombination for use in generating *C. elegans* transgenic animals by bombardment. *Plasmid, 62*, 140-145.

# **Buffer Preparation**

✓ <u>Before starting</u>: Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **DNA** Wash Buffer concentrate.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

## **Protocol**

Add 2-7 volumes of **DNA Binding Buffer** to each volume of DNA sample (see table below). Mix briefly by gently inverting the tube.

Application	DNA Binding Buffer : Sample	Example
Plasmid, genomic DNA (>2 kb)	2:1	200 µl : 100 µl
PCR product, DNA fragment	5 : 1	500 μl : 100 μl
ssDNA (e.g. cDNA, M13 phage1)	7 : 1	700 μΙ : 100 μΙ

Use any of the following two procedures to process samples.

## Centrifuge <sup>2</sup>

- 1. Place the **Zymo-Spin<sup>™</sup> VI Column** inside a 50 ml conical tube. Transfer the prepared sample mixture (see above) into the **Zymo-Spin<sup>™</sup> VI Column**<sup>3</sup>.
- 2. Centrifuge for 5 minutes at 3,000 x  $g^2$ . Discard the flow-through.
- 3. Add 10 ml DNA Wash Buffer to the Zymo-Spin™ VI Column. Centrifuge the Zymo-Spin™ VI Column for 5 minutes.
- 4. Transfer the **Zymo-Spin™ VI Column** into a new 50 ml conical tube. Add 2-3 ml **DNA Elution Buffer**<sup>4</sup> or water<sup>5</sup> directly to the column matrix. Wait for one minute to ensure that the column matrix has been fully hydrated. Centrifuge for 3 minutes to elute the DNA.

Ultra-pure DNA is now ready for use.

#### Vacuum

- Attach Zymo-Spin™ VI column to a suitable vacuum manifold (see image on page 2). Transfer the prepared sample mixture (see above) into the Zymo-Spin™ VI Column³.
- 2. Turn on vacuum source and wait for sample to clear from the column.
- 3. Add 5 ml **DNA Wash Buffer** to the **Zymo-Spin™ VI Column**. Repeat wash step. After washing, leave the vacuum source "on" for an additional 5 minutes to remove any residual **Wash Buffer** from the column.
- 4. Transfer the **Zymo-Spin™ VI Column** into a new 50 ml conical tube. Add 2-3 ml **DNA Elution Buffer**<sup>4</sup> or water<sup>5</sup> directly to the column matrix. Wait for one minute to ensure that the column matrix has been fully hydrated. Centrifuge at 3,000 x *g* for 3 minutes to elute the DNA.

Ultra-pure DNA is now ready for use.

#### Notes:

- <sup>1</sup> Centrifuge phage-infected bacterial culture at 8,000 x g for 1 minute prior to mixing an aliquot of the phage-containing supernatant with the **DNA Binding Buffer**.
- <sup>2</sup> All centrifugation steps should be performed at 3,000 x *g*.
- The sample capacity of the column is 15 ml. Therefore, it may be necessary to load a column multiple times if a sample has a volume larger than 15 ml.
- <sup>4</sup> **DNA Elution Buffer**: 10 mM Tris-HCl, pH 8.5, 0.1 mM EDTA
- <sup>5</sup> Elution of DNA from the column is dependent on pH and temperature. If water is used, make sure the pH is >6.0. Waiting an additional minute prior to elution may improve the yield of larger (> 6 kb) DNA. For even larger DNA (> 10 kb), the total yield may be improved by eluting the DNA with 60-70 °C DNA Elution Buffer.

## **Troubleshooting**

### Low Recovery

### • Improperly Prepared/Stored DNA Wash Buffer

Make sure ethanol has been added to the **DNA Wash Buffer** concentrate. Cap the bottle tightly to prevent evaporation over time.

#### Addition of DNA Elution Buffer

Add elution buffer directly to the column matrix and not to the walls of the column. Elution buffer requires contact with the matrix for at least 1 minute for large DNA  $\geq$  10 kb.

#### • Incomplete Elution

- DNA elution is dependent on pH, temperature, and time. For large genomic DNA (≥ 50 kb), apply heated elution buffer (60-70 °C) and incubate for several minutes prior to elution.
- 2. Sequential elutions may be performed for quantitatively higher recovery but lower final DNA concentration. This is recommended for DNA ≥ 10 kb.

# Low A<sub>260</sub>/A<sub>230</sub> Ratios

#### Column Tip Contaminated

When removing the column from the collection tube, be careful that the tip of the column does not come into contact with the flowthrough. Trace amounts of salt from the flowthrough can contaminate a sample resulting in low  $A_{260}/A_{230}$  ratios. Ethanol contamination from the flowthrough can also interfere with DNA elution. Zymo-Spin<sup>TM</sup> columns are designed for complete elution with no buffer retention or carryover.

# Following Clean-up with the DCC™, Multiple Bands Appear in an Agarose Gel

#### Acidification of DNA Loading Dye

Most loading dyes do not contain EDTA and will acidify (pH  $\leq$  4) over time due to some microbial growth. This low pH is enough to cause DNA degradation. Therefore, if water is used to elute the DNA, 6X Loading Dye containing 1 mM EDTA is recommended.

# **Ordering Information**

Product Description	Catalog No.	Kit Size (Preps.)
DNA Clean & Concentrator <sup>TM</sup> -5 (for purification of up to 5 µg DNA per prep.) Supplied with uncapped columns	D4003 D4004	50 200
DNA Clean & Concentrator <sup>TM</sup> -5 (for purification of up to 5 µg DNA per prep.) Supplied with capped columns	D4013 D4014	50 200
ZR-96 DNA Clean & Concentrator TM-5 (for 96-well purification of up to 5 μg DNA per well)	D4023 D4024	2 x 96 4 x 96
DNA Clean & Concentrator <sup>TM</sup> -25 (for purification of up to 25 μg DNA per prep.) Supplied with uncapped columns	D4005 D4006	50 200
DNA Clean & Concentrator <sup>TM</sup> -25 (for purification of up to 25 μg DNA per prep.) Supplied with capped columns	D4033 D4034	50 200
DNA Clean & Concentrator <sup>TM</sup> -100 (for purification of up to 100 μg DNA per prep.)	D4029 D4030	25 50
DNA Clean & Concentrator <sup>TM</sup> -500 (for purification of up to 500 μg DNA per prep.)	D4031 D4032	10 20
Oligo Clean & Concentrator <sup>TM</sup> (for purification of up to 5 μg of oligonucleotides per prep.)	D4060 D4061	50 200
Genomic DNA Clean & Concentrator <sup>TM</sup> (for purification of up to 10 μg genomic DNA per prep.)	D4010 D4011	25 100

# Refer to Page 3 for column design specifics in each kit.

For Individual Sale	Catalog No.	Size
DNA Binding Buffer	D4003-1-L D4004-1-L	50 ml 100 ml
DNA Wash Buffer (concentrate)	D4003-2-24 D4003-2-48	24 ml 48 ml
DNA Elution Buffer	D3004-4-16 D3004-4-50	16 ml 50 ml
Zymo-Spin™ VI Columns	C1013-10 C1013-20	10 columns 20 columns

# Popular Products From Zymo Research

Product	Description	Kit Size (Preps.)	Catalog No. (Format)
	DNA Clean-up, Concentration & Recovery		
DNA Clean & Concentrator™-5	Clean and concentrate up to 5 µg DNA into ≥6 µl elution volume in as little as 2 minutes with no wash residue carryover.	50 200 50 200	D4003 (uncapped) D4004 (uncapped) D4013 (capped) D4014 (capped)
DNA Clean & Concentrator™-25	Clean and concentrate 25 µg of DNA into ≥25 µl elution volume in as little as 2 minutes with no wash residue carryover.	50 200 50 200	D4005 (uncapped) D4006 (uncapped) D4033 (capped) D4034 (capped)
ZR-96 DNA Clean &	Quick (30 minute), high throughput recovery of up to 5 µg pure DNA into 10-15 µl minimum elution volume allows for highly concentrated DNA.	2 x 96	D4023
Concentrator™-5		4 x 96	D4024
Genomic DNA Clean &	Quick (5 minute) clean-up of up to 10 µg high molecular weight DNA (≥ 20kb - 200 kb) from any enzymatic reaction or impure preparation without precipitations.	25	D4010
Concentrator™		100	D4011
Zymoclean™ Gel DNA Recovery Kit	Purify DNA from high and low-melting agarose gels in minutes.	50 200 50 200	D4001 (uncapped) D4002 (uncapped) D4007 (capped) D4008 (capped)
ZR-96 Zymoclean™ Gel	High-throughput DNA purification from high and low-melting agarose gels.	2 x 96	D4021
DNA Recovery Kit		4 x 96	D4022
Zymoclean™ Large Fragment	Purify high molecular weight DNA (≥ 20 kb - 200 kb) from high and low-melting agarose gels in minutes.	25	D4045
DNA Recovery Kit		100	D4046
OneStep™ PCR	Fast, one step procedure for removal of PCR inhibitors such as polyphenolics, humic/fulvic acids, melanin, etc. for successful PCR and other downstream applications.	50	D6030
Inhibitor Removal Kit		2 x 96	D6035

Plasmid DNA Purification			
Zyppy™ Plasmid Miniprep Kit	Pellet-Free™ plasmid DNA purification in less than 10 minutes. Recover up to 25 μg DNA in as low as 30 μl.	50 100 400	D4036 D4019 D4020
Zyppy™-96 Plasmid Miniprep	The fastest and simplest high-throughput method for plasmid purification. Magnetic bead format available for automated liquid handling platforms.	2 x 96 4 x 96 8 x 96 2 x 96 4 x 96 8 x 96	D4041 (spin plate) D4042 (spin plate) D4043 (spin plate) D4100 (magnetic bead) D4101 (magnetic bead) D4102 (magnetic bead)
Zyppy™ Plasmid Midiprep Kit	Pellet-Free™ plasmid DNA purification in 15 minutes in a 150 µl minimum elution volume.	25 50	D4025 D4026
ZR Plasmid MiniPrep™- <i>Classic</i>	Plasmid DNA purification in minutes: (alkaline lysis/spin column format for low 30 µl elution volume).	100 400 800	D4015 D4016 D4054

Genomic DNA Purification			
<i>Quick-gDNA</i> ™ MiniPrep	Easy purification from whole blood, plasma, serum, body fluids, buffy coat, tissue, swabs or cultured cells ≥15 minutes <u>without</u> the use of Proteinase K or organic denaturants.	50/200 50/200	D3006/D3007 uncapped) D3024/D3025 (capped)
ZR Genomic DNA™-Tissue MiniPrep	igh quality DNA purification from solid tissues (e.g., tail snips, ear punches, adipose tissue, etc.), body fluids, cultured cells, buccal cells, FFPE tissues, hair, and other biological sources using Proteinase K and Fast.	50/200 50/200	D3050 D3051
Environmental DNA Purification Kits	Unique BashingBead™ technology allows isolation of DNA from samples refractory to conventional lysis procedures including tough-to-lyse tissues, soil samples, feces, plants, seeds, insects, bacteria, yeast, filamentous fungi, unicellular and filamentous algae, and protozoa	Spin Column & 96-well Plate	Visit website for a comprehensive list

RNA Purification			
RNA Clean & Concentrator™-5	Clean and concentrate up to 5 µg RNA into ≥6 µl elution volume in as little as 5 minutes with no wash residue carryover.	50 200	R1015 R1016
Direct-Zol™ RNA MiniPrep	Quick, spin column purification of high-quality (DNA-free) total RNA <i>directly</i> from <i>TRI-Reagent</i> ® or similar acid-guanidinium-phenol based reagents (TRIzol®, RNAzol®, QIAzol®, TriPure, RNA-Bee <i>etc.</i> ).	50 200	R2051 R2053
ZR RNA MiniPrep	Quick (15 minute) RNA isolation from a variety of sources using Fast-Spin column technology without the use of organic denaturants	50 200	R1064 R1065

Epigenetics Products From Zymo Research



Product	Description	Kit Size	Cat. No. (Format)
	Bisulfite Kits for DNA Methylation Detection		
EZ DNA Methylation™ Kit	For the conversion of unmethylated cytosines in DNA to uracil via the <a href="https://enaturation.com/chemostration">chemostration</a> of DNA and a specially designed CT Conversion Reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5001/D5002 (column) D5003 (shallow-well plate) D5004 (deep-well plate) D5040 (magnetic bead)
EZ DNA Methylation- Gold™ Kit	For the fast (3 hr.) conversion of unmethylated cytosines in DNA to uracil via <a href="heat/chemical-denaturation">heat/chemical-denaturation</a> of DNA and a specially designed CT Conversion Reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5005/D5006 (column) D5007 (shallow-well plate) D5008 (deep-well plate) D5042 (magnetic bead)
EZ DNA Methylation- Direct™ Kit	Simple and reliable DNA bisulfite conversion directly from blood, tissue (FFPE/LCM), and cells without the prerequisite for DNA purification in as little as 4-6 hrs. The increased sensitivity of this kit makes it possible to amplify bisulfite converted DNA from as few as 10 cells or 50 pg DNA. Magnetic bead format for adaptation to automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5020/D5021 (column) D5022 (shallow-well plate) D5023 (deep-well plate) D5044 (magnetic bead)
EZ DNA Methylation- Lightning™ Kit	Complete bisulfite conversion in about an hour using a unique liquid format conversion reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5030/D5031 (column) D5032 (shallow-well plate) D5033 (deep-well plate) D5046 (magnetic bead)
EZ DNA Methylation- Startup™ Kit	Designed for the first time user requiring a consolidated product to perform DNA methylation analysis. Includes technologies for sample processing, bisulfite treatment of DNA, and PCR amplification of "converted" DNA for methylation analysis.	1 Kit	D5024
	Methylated DNA Standards		
Universal Methylated Human DNA Standard	Human (male) genomic DNA having all CpG sites methylated. To be used for the evaluation of bisulfite-mediated conversion of DNA. Supplied with a control primer set.	1 Set	D5011
Universal Methylated Mouse DNA Standard	Mouse (male) DNA having all CpG sites methylated. To be used for the evaluation of bisulfite-mediated conversion of DNA. Supplied with a control primer set.	1 Set	D5012
	Region-Specific DNA Methylation Screening		
<i>OneStep</i> qMethyl™ Kit	Single step real-time PCR procedure for bisulfite-free determination of DNA methylation status. Available without fluorescent dye for probe-based detection (Lite).	1 x 96 Rxns. 1 x 96 Rxns.	D5310 D5311 (Lite)
<i>OneStep</i> qMethyl™ Array	Premade 96-well assay for bisulfite-free determination of region-specific DNA methylation assessment in the promoter region of any one of the following prominent tumor suppressor genes: RASSF1, RARB, CDKN2A (p16), MGMT, or CCND2.	1 x 96 Rxns.	D5312
	Epigenetics Services		
comprehensive DNA methylat Services for Hydroxymethyl	with our 5-mC Analysis platforms that combine Zymo's well-established bisulfite technologic ion analysis services available. ated DNA Analysis	es with next-genera	ntion sequencing for the most
Novel genome-wide 5-hmC au sensitivity of 5-hmC detection	nalysis platform featuring cutting-edge 5-hmC DNA enrichment, library prep, and next-gene in genome-wide context.	ration sequencing	echnologies to ensure the
	Hydroxymethylation Detection		
Quest 5-hmC™ DNA Enrichment Kit	Featuring J-base binding protein (JBP) for the specific enrichment of 5-hmC containing DNA, the consolidated workflow makes the procedure reliable for robust analysis of multiple samples.	25 Rxns. 50 Rxns.	D5420 D5421
Quest 5-hmC™ DNA ELISA Kit	Streamlined workflow for both the direct and relative quantitation of 5-hmC, in a global genomic context, with a robust colorimetric readout.	1 x 96 Rxns. 2 x 96 Rxns.	D5425 D5426
Anti-5- Hydroxymethylcytosine Polyclonal Antibody	Polyclonal antibody has been engineered to maximize sensitivity to low amounts of hydroxymethylated gDNA while minimizing crossreactivity with unmodified or methylated cytosine residues. The antibody is suitable for use in ELISA, IP, and immunohistochemical labeling.	50 μg 200 μg	A4001-50 A4001-200
DNA Degradase™ DNA Degradase Plus™	Whole genomic DNA can be treated with these enzyme cocktails for processing to individual nucleotides (Degradase™) or nucleosides (Degradase Plus™) for interrogation in chromatographic and spectroscopic methods including TLC, LC/MS, MALDI-TOF, and more.	500 U 2000 U 250 U 1000 U	E2016 E2017 E2020 E2021
	Other		
Zymo <i>Taq</i> ™ DNA Polymerase	Zymo Taq <sup>™</sup> "hot start" DNA Polymerase is specifically designed for the amplification of "difficult" DNA templates including: bisulfite-treated DNA for methylation detection. The product generates specific amplicons with little or no by-product formation.	50 Rxns. 200 Rxns	E2001/E2001 (system) E2003/E2004 (premix)
Methylated-DNA IP Kit	IP with a highly specific anti-5-methylcytosine monoclonal antibody. Designed for the enrichment of 5-methylcytosine-containing DNA from any pool of fragmented genomic	10 Rxns.	D5101

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com