

INSTRUCTION MANUAL

RNA Clean & Concentrator[™]-100

Catalog Nos. R1019

Highlights

- Quick, 15 minute recovery of ultra pure RNA (≥17 nt) from enzymatic reactions, aqueous phase (following Trizol® extraction) and other sources.
- High-quality RNA eluted in ≥100 µl is ready for reverse transcription, microarray, sequencing, etc.

Contents

Product Contents	
Specifications	
Product Description	2
Buffer Preparation	
Protocol	
Appendices	
Ordering Information	
Related Products	

For Research Use Only Ver. 1.0.3

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product please call 1-888-882-9682.

Notes:

¹ Before use, add 24 ml 100% ethanol (26 ml 95% ethanol) to the 6 ml **RNA Wash Buffer** concentrate (R1019).

Product Contents

RNA Clean & Concentrator [™] -100 (Kit Size)	R1019 (25 Preps.)
RNA Binding Buffer	100 ml
RNA Prep Buffer	10 ml
RNA Wash Buffer¹ (concentrate)	6 ml
DNase/RNase-Free Water	10 ml
Zymo-Spin [™] V-E Columns w/ Reservoir	25
Collection Tubes	50
Instruction Manual	1

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Storage Temperature – Store all kit components (*i.e.*, buffers, columns) at room temperature.

Specifications

- **Sample Sources** DNase I treated RNA, *in vitro* transcription products, the aqueous phase following TRIzol®/chloroform or similar² extraction (page 4).
- RNA Size Limits From 17 nt to ~23 kb.
- **RNA Purity** High quality RNA ($A_{260}/A_{280} > 1.8$, $A_{260}/A_{230} > 1.8$) suitable for reverse transcription, microarray, sequencing etc.
- RNA Recovery Up to 250 µg RNA can be eluted into as little as ≥100 µl RNase-free water allowing for a highly concentrated sample.
- RNA Storage RNA is eluted with RNase-free water and can be stored at ≤-70 °C. The addition of RNase inhibitors is optional but highly recommended for prolonged storage.
- **Equipment Needed** Vacuum manifold, microcentrifuge.

TRI Reagent®, RNAzol®, QIAzol®, TriPure™, TriSure™ and other acid-guanidinium-phenol reagents.

² Compatible with: TRIzol[®],

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

TRI Reagent[®], TRIzol[®] and RNAzol[®] (Molecular Research Center, Inc.), QIÁzol[®] (Qiagen GmbH), TriPure[™] (Roche Diagnostics Operations, Inc.), TriSure[™] (Bioline Ltd.), RNA*latel*[®] (Ambion, Inc.).


Product Description


RNA Clean & Concentrator[™]-100 provides a simple and reliable method for the rapid preparation of up to ~250 µg of high-quality RT-PCR-ready RNA. This simple procedure is based on the use of a unique single-buffer system and Clean-Spin[™] column technology that allows for selective recovery of total RNA (> 17 nt), large RNAs (> 200 nt), and/or small RNAs (17-200 nt).

The procedure is easy: Add binding buffer and ethanol to your sample, then bind, wash and elute ultra pure RNA. The RNA can be eluted from the **Zymo-Spin**TM **V-E Column** in as little as \geq 100 μ I of RNase-free water. The highly-concentrated, purified RNA is suitable for all subsequent analyses and molecular manipulations.

The entire procedure typically takes about 15 minutes.

For Assistance, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.





Concentration of diluted RNA samples. RNA was eluted with 20 µl RNase-free water (n = 3, total input = 1 µg RNA). Purification of small and large RNAs into separate fractions. RNA Clean & Concentrator™ allows for purification of total RNA (> 17 nt), large RNAs (> 200 nt), and/or small RNAs (17-200 nt).

Note:

For purification of DNA see the **DNA Clean & Concentrator**™-**5** and **-25** (Catalog Nos. D4013, D4014, D4033, D4034).

Make sure guidelines are followed to ensure the RNA isolation procedure is performed in an RNase-free environment.

Buffer Preparation

Before starting, add 48 ml 100% ethanol (52 ml 95% ethanol) to the 12 ml **RNA Wash Buffer** concentrate (R1017) or 96 ml 100% ethanol (104 ml of 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate (R1018).

Protocol

All centrifugation steps should be performed at $10,000 - 16,000 \times g$. RNA species ≥ 17 nt will be recovered.

1. Add 2 volumes **RNA Binding Buffer** to each sample¹ and mix.

Example: Mix 1 ml buffer and 0.5 ml sample.

2. Add an equal volume of ethanol (95-100%) and mix.

Example: Add 1.5 ml ethanol.

- 3. Transfer the sample to the **Zymo-Spin**[™] **V-E Column w/ Reservoir** mounted on a vacuum manifold and start vacuum^{1,2}. Discard the flow-through.
- 4. Add 400 μl **RNA Prep Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
- 5. Add 400 μl **RNA Wash Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
- Add 400 µl RNA Wash Buffer to the column and centrifuge for 2 minutes to ensure complete removal of the wash buffer. Transfer the column carefully into an RNasefree tube (not provided).
- 7. Add ≥100 µl **DNase/RNase-Free Water** directly to the column matrix and centrifuge for 30 seconds.

The eluted RNA can be used immediately or stored at -70°C.

Notes:

- ¹ Set vacuum source at ≥500 mm Hg.
- ² At this point, samples can be in-column DNase treated (page 4).

DNase I treatment

The DNase digestion procedure can be performed using the **DNase I Set** (E1010)¹.

1. For each sample to be treated, prepare **DNase I reaction mix** in an RNase-free tube (not provided). Mix well by gentle inversion:

RNA sample (≤250 µg)

volume adjusted with water or TE buffer 160 µl 20 µl DNase I 20 µl DNA Digestion Buffer 20 µl 200 µl

2. Incubate at room temperature (20-30°C) for 15 minutes. Then start with RNA purification (page 3, step 1).

In-Column DNase I treatment

All centrifugation steps should be performed at $10,000 - 16,000 \times g$.

- 1. Following the RNA binding step (page 3, step 3), prewash the column with 400 µl **RNA Wash Buffer**. Centrifuge for 30 seconds. Discard the flow through.
- 2. For each sample to be treated, prepare **DNase I reaction mix** in an RNase-free tube (not provided). Mix well by gentle inversion:

DNase I 20 μl **DNA Digestion Buffer** 300 μl

3. Add 320 µl **DNase I reaction mix** directly to the column matrix. Incubate at room temperature (20-30°C) for 15 minutes. Then centrifuge the column for 30 seconds and discard the flow-through. Continue with RNA purification (page 3, step 4).

RNA purification from aqueous phase after TRIzol® extraction

- 1. Following Trizol®/chloroform or similar² extraction, carefully transfer the upper aqueous phase into an RNase-free tube (not provided).
- 2. For each volume of the aqueous phase (as measured or estimated), add 1 volume ethanol (95-100%) and mix.
- 3. Then continue with purification (page 3, step 3).

Notes:

¹ Prior to use, reconstitute the lyophilized **DNase I** as indicated on the vial. Store frozen aliquots.

Unit definition - one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A₂₆₀ units/min/ml of reaction mixture at 25°C.


² Compatible with: TRIzol®, TRI Reagent®, RNAzol®, QIAzol®, TriPure™, TriSure™ and other acid-guanidinium-phenol reagents.

Ordering Information

Product Description	Catalog No.	Kit Size
RNA Clean & Concentrator™-100	R1019	25 Preps.
RNA Clean & Concentrator™-25	R1017 R1018	50 Preps. 100 Preps.
RNA Clean & Concentrator™-5	R1015 R1016	50 Preps. 200 Preps.
RNA Clean & Concentrator [™] -5 with DNase I Set	R1013 R1014	50 Preps. 200 Preps.

For Individual Sale	Catalog No.	Amount
RNA Binding Buffer	R1013-2-25 R1013-2-50 R1013-2-100 R1013-2-1000	25 ml 50 ml 100 ml 1000 ml
RNA Prep Buffer	R1060-2-10 R1060-2-25 R1060-2-100	10 ml 25 ml 100 ml
RNA Wash Buffer (concentrate)	R1003-3-6 R1003-3-12 R1003-3-24 R1003-3-48	6 ml 12 ml 24 ml 48 ml
Zymo-Spin V-E Columns w/ Reservoir	C1024-25	25
Collection Tubes	C1001-50 C1001-500 C1001-1000	50 500 1000
DNase/RNase-Free Water	W1001-1 W1001-4 W1001-6 W1001-10	1 ml 4 ml 6 ml 10 ml

DNA PURIFICATION




Purify DNA from PCR & other sources

DNA Clean & Concentrator™ (DCC™)

- ✓ Recovery of ultra-pure DNA that is free of salts and contaminants.
- ✓ Small (≥6 µI) elution volume.
- ✓ DNA is ideal for ligation, PCR, Next-Gen sequencing, etc.

Product	Size (Cat. No.)
DNA Clean & Concentrator™-5	50 Preps. (D4013) 200 Preps. (D4014)
ZR-96 DNA Clean & Concentrator™-5	2 x 96 Preps. (D4023) 4 x 96 Preps. (D4024)
Genomic DNA Clean & Concentrator™	25 Preps. (D4010) 100 Preps. (D4011)




Boost DNA recoveries from agarose gels to >80%

Zymoclean™ Gel DNA Recovery

- ✓ Rapid (15 min.) recovery of ultra-pure DNA from agarose gels in ≥6 μl.
- ✓ Ultra-pure DNA ideal for DNA ligation, sequencing, etc.
- √ Format also available for large DNA >20 kb.

Product	Size (Cat. No.)
Zymoclean™ Gel DNA Recovery Kit	50 Preps. (D4001) 200 Preps. (D4002)
Zymoclean™ Large Fragment DNA Recovery Kit	25 Preps. (D4045) 100 Preps. (D4046)



ladder DNA fragments.

Recover transfection-quality plasmid DNA directly from culture

Zyppy™ Plasmid Prep Kits

- ✓ The fastest, simplest method available for purifying high quality plasmid DNA from E. coli.
- ✓ Pellet-Free[™] procedure omits conventional cell-pelleting and resuspension steps.
- √ Transfection quality plasmid DNA directly from culture in under 15 minutes.



Product	Size (Cat. No.)
Zyppy [™] Plasmid Miniprep Kit	50 Preps. (D4036) 100 Preps. (D4019) 400 Preps. (D4020) 800 Preps. (D4037)



What is Clean-Spin[™] Technology? The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Get RNA <u>directly</u> from TRIzol® without phase separation


Direct-zol™ RNA

- ✓ For purification of high-quality small and large RNA directly from TRIzol®, TRI Reagent®, or similar.
- ✓ Bypasses phase separation and precipitation procedures allowing for unbiased recovery of miRNA

Product	Size (Cat. No.)
Direct-zol™ RNA MiniPrep	50 Preps. (R2050) 50 Preps. (R2051)* 200 Preps. (R2052) 200 Preps. (R2053)*
96-well and MagBead formats also available!	

DNase Lincluded in all kits

Supplied with TRI-Reagent®




Isolate DNA-free RNA from 1 to 107 cells in minutes

Quick-RNA™

- ✓ Isolation of total, large, or small RNA You decide!
- ✓ Ultra clean, high-quality RNA from a single cell to 10⁷ cells.
- DNA-free RNA ideal for any downstream application DNase I included.

Product	Size (Cat. No.)
Quick-RNA™ MicroPrep	50 Preps. (R1050) 200 Preps. (R1051)
Quick-RNA™ MiniPrep	50 Preps. (R1054) 200 Preps. (R1055)
ZR-96 Quick-RNA™	2 x 96 Preps. (R1052) 4 x 96 Preps. (R1053)



Purify RNA from enzymatic and labeling reactions in 5 minutes

RNA Clean & Concentrator™

- ✓ Recover ultra-pure RNA in small (≥6 μl) elution volumes.
- ✓ Compatible with TRIzol®, phenol, choloform, and RNase inhibitors (RNAlater®).
- ✓ RNA is ideal for RT-PCR, q-PCR, hybridization, arrays, RNA interference, etc.

Product	Size (Cat. No.)
RNA Clean & Concentrator™-5	50 Preps. (R1015) 200 Preps. (R1016)
RNA Clean & Concentrator™-25	50 Preps. (R1017) 100 Preps. (R1018)
ZR-96 RNA Clean & Concentrator™	2x96 well plates (R1080)
DNA-Free RNA Kit™	50 Preps. (R1013) 200 Preps. (R1014)



The following are trademarks of other companies: pGEM®, Promega Corp.; TRIzol® and TRI Reagent®, Molecular Research Center, Inc.; DH5® and DH10B™, Life Technologies, Inc.



The Beauty of Science is to Make Things Simple