

Flow Cytometric Screening of a Kinase Inhibitor Library for Apoptosis Induction

Julie M. Rumble and David L. Hoffman


Bioanalysis & Assay Development Services Group, Cayman Chemical, Ann Arbor, MI

KEY FINDING Multiparametric screening of kinase inhibitor cytotoxicity distinguishes multiple modes of cell death.

Introduction

Traditional methods of detecting cytotoxicity of potential therapeutic compounds have relied on less sensitive methods, which are amenable to high throughput. These include detecting LDH release, tetrazolium salt reduction to formazan, and ATP. These methods can be useful, but finer understanding of the mechanism of cell death may be more helpful in developing effective therapeutics.

Fluorescent probes can detect many different components of cell death and have the advantage of being able to be multiplexed, given the proper analysis platform. Annexin V binds to externalized phosphatidylserine, a marker of apoptosis. DAPI is a nucleic acid dye, which is only internalized after the plasma membrane is compromised in necrotic cells. TO-PRO[®]-3 is trafficked into the cell early in apoptosis through pannexin channels. TMRE is fluorescent only in healthy mitochondria.




Flow cytometry has not traditionally been used as a screening tool due to low throughput. That is, samples are analyzed individually, and automation is frequently unreliable, thus demanding high user input. The MACSQuant[®] X has been developed to counter these issues and make flow cytometry significantly higher throughput. This flow cytometer can analyze as little as 5 µl of sample reliably, can automatically sample 384-well plates, and maintains the multiplexability that makes flow cytometry an attractive tool for obtaining maximal information from each sample.

Methods

Jurkat cells were plated at 4×10^5 /well in two 96-well plates in complete medium and treated overnight with 1 µM of compound from the Kinase Screening Library (Cayman Item No. 10505). Approximately 160 compounds are in this library, and the remainder of the wells were used as controls. DMSO was used as a vehicle control, and staurosporine (Cayman Item No. 81590) was used as a positive control. Cells were stained as described in the kit booklet for the Early Apoptosis Detection Assay Kit (Cayman Item No. 601360), using annexin V FITC, TMRE, and TO-PRO[®]-3. After staining, cells were washed, resuspended in PBS with DAPI, and transferred to 384-well plates in quadruplicate. Cells were analyzed in the 384-well plates using the MACSQuant[®] X flow cytometer with the following settings: Collect 15 µl sample, high speed, fast wash, shake the plate every 6 wells; collect data in channels V1, B1, B2, and R1. Using these settings, each 384-well plate was sampled in about 2.5 hours.

Cell Staining and MACSQuant[®] X Flow Cytometry Workflow



Results

Analysis of Flow Cytometric Data




Figure 1. A. An initial gate was drawn to enumerate necrotic cells (DAPI+). B. DAPI-negative cells were further gated on annexin V-positive and -negative cells. C. The annexin V-positive population was then gated using forward scatter to differentiate cells from apoptotic bodies and debris. D. The annexin V-negative cells were visualized in a histogram of TO-PRO[®]-3, in which positive cells were described as early apoptotic and negative cells were live.

Heat Map Reveals Potentially Cytotoxic Compounds

Item No.	Item Name	Live	Item No.	Item Name	Live	Item No.	Item Name	Live
84.88	TG003	85.13	13341	LY364947	85.68	10010236	ML-9	85.76
11022	BHFF 1120	85.25	13812	17β-hydroxy Wortmannin	85.73	10010248	AG-17	85.76
11793	SB-505124	86.25	9000988	PD 160326	78.85	10010312	AG-82	85.38
85.78	CHIR99021	79.43	10793	AG-879	84.23	10010466	SP 600125	85.05
85.23	NSC 663284	115.69	11569	GSK1059615	84.23	10011249	CAY10576	85.00
13325	Iso-Olomoucine	84.28	13067	SB 203580	84.23	10008005	OSU03012	80.93
13601	ZM447439	81.23	13218	Bisindolylmaleimide I	85.23	10007707	AS-605240	84.98
62575	N,N-Dimethylsphingosine (d18:1)	82.90	13318	KN-62	85.78	10009222	Sphingosine Kinase Inhibitor 2	85.20
10459	PKC 412	81.05	13344	SB 203580 (hydrochloride)	84.60	10010237	Triciribine	85.18
11029	SMI-4a	84.88	13388	CAY10626	84.10	10010249	H-8 (hydrochloride)	86.58
13811	INK128	79.85	10004914	O-1518	79.88	10010313	AG-99	86.30
13123	RO	85.88	11444	L-NA-PP1	85.40	10010541	Leuco-Sphingosine (d18:1)	85.73
13305	D 4476	85.65	11609	Ruxolitinib	84.20	Control	Control	85.15
13332	(S)-Glycyl-H-1152 (hydrochloride)	84.60	13108	VX-702	83.83	10010301	CAY10554	82.50
13622	AG-041164	84.43	13299	Bisindolylmaleimide IV	84.93	10007907	Sphingosine (d18:1)	85.00
70920	LY294002	82.88	13319	KN-93	84.38	10009366	Picetannolol	85.18
10460	Doramapimod	83.95	13371	CAY10621	83.93	10010238	Erbstatin analog	85.45
13140	CAY10657	81.63	13873	SU6668	83.68	10010265	LFM-A13	86.25
12076	Canertinib (hydrochloride)	84.45	10005583	V-27632 (hydrochloride)	82.95	10010314	AG-213	85.90
13139	Imatinib (mesylate)	84.83	10997	Torin 1	81.25	10010556	H-89 (hydrochloride)	85.83
13308	NU 7076	84.75	11658	Necrostatin-1	84.00	10011251	TWS119	81.93
13333	Bisindolylmaleimide VIII (acetate)	86.68	13109	Emodin	84.63	Control	Control	85.15
13641	NVP-AEW541 (hydrochloride)	84.38	13300	Bisindolylmaleimide V	85.23	10008131	JNJ-10198409	55.63
70970	U-0126	88.58	13322	CGP 57380	83.25	10009557	SC-1	75.60
10461	Paclitaxel	59.48	13576	HA-201636	80.28	10010239	Kemupallone	86.68
Control	Control	84.85	18218	PHA-767491	84.28	10010267	SC-514	86.60
13031	SB-431542	84.98	10006148	Leleamine	82.95	10010315	AG-183	86.35
13159	Sunitinib (maleate)	83.50	10006726	PD 98059	84.63	10010559	HA-1077 (hydrochloride)	85.98
13311	Gö 6983	85.30	10009052	AS-252424	85.60	10011255	NSC 210902	85.25
13334	Bisindolylmaleimide IX (mesylate)	86.80	10010175	AS-604850	85.10	Control	Control	80.50
13643	PP242	81.60	10010244	AG-1478	85.28	10008614	Leleamine (hydrochloride)	85.10
81590	Staurosporine	0.90	10010309	RG-13022	84.98	10009569	(R)-Roscovitine	85.60
10483	Erlotinib	83.80	10010399	SB 202190	83.95	10010240	Olomoucine	86.20
11314	Chloroethylamine (chloride)	85.38	10011246	WH-P131	82.58	10010275	Ajagennin	86.15
13032	PD 173074	85.53	10012591	CCT018159	79.98	10010329	Levendustin C	85.95
13166	Gefitinib	84.98	10006727	PD 169316	84.85	10010568	AG-370	85.63
13312	H-9 (hydrochloride)	85.15	10009078	CAY10505	84.73	10011256	CAY10577	80.03
13337	ST636	85.30	10010177	PI3-Kinase α Inhibitor 2	85.25	Control	Control	80.20
13653	ABT-689	84.45	10010246	SB-126763	86.68	10010242	AG-494	84.75
Control	Control	84.85	10010310	RG-14620	85.45	10009644	Sorafenib	85.40
10565	NVP-BE2235	74.73	10010400	CAY10571	85.30	10010242	AG-494	85.80
11445	Tunicamycin	67.13	10011247	CAY10574	82.03	10010300	AG-18	86.20
13033	Valproic Acid (sodium salt)	83.83	10012600	Myricetin	78.50	Control	Control	81.03
13198	PP2	84.98	10007349	TGK-221	85.10	10010591	Wortmannin	85.45
13314	Indirubin-3'-monoxime	84.68	10009209	PI-103	85.45	10011264	CAY10578	83.60
13388	SU6656	84.85	10010233	CAY10567	85.43	Control	Control	79.20
13687	CAY10622	84.33	10010247	SB-415286	86.93	10008618	Lauric Acid Leleamine	84.45
9009980	AS-605240 (potassium salt)	82.98	10010312	AG-225	85.28	10010043	CAY10561	85.05
10735	Phthalazinone pyrazole	83.88	10010422	Nilotinib	85.33	10010243	AG-825	85.10
11491	AZD 7762	35.73	Control	Control	81.60	10010302	DRB	84.70
13034	PD 0325901	84.08	10527	Necrostatin-5	78.90	10010375	S-lodotubercidin	84.78
13242	3-Methyladenine	85.50	10007653	(S)-1152 (hydrochloride)	83.10	10010592	AG-1296	82.83
13317	NU 6102	84.08	10009210	PIK-75 (hydrochloride)	2.29	10012431	PD 184161	82.73
			Control	Control	80.93	Control	Control	80.93

Figure 2. The percentage of the total population falling into the "live" (DAPI-, annexin V-, TO-PRO[®]-3-) category was used to flag cytotoxic compounds (blue; live, red; dead). The average of the DMSO and staurosporine-treated controls was 83% and 1.6%, respectively.

Further Analysis of Potentially Cytotoxic Compounds

Item No.	Item Name	Live	Early Apoptotic	Necrotic	Apoptotic
10461	Paclitaxel	59.48	15.075	10.975	3.91
81590	Staurosporine	0.90	12.425	68.575	3.47
10565	NVP-BE2235	74.73	8.295	11.1	1.08
11445	Tunicamycin	67.13	8.47	13.825	4.01
11491	AZD 7762	35.73	18.625	26.575	4.96
13812	17β-hydroxy Wortmannin	57.13	7.3525	29.975	0.69
11569	GSK1059615	76.18	6.4875	11.825	1.55
10009210	PIK-75 (hydrochloride)	2.29	20.425	55.775	10.83
10010248	AG-17	76.13	3.4775	16.4	0.47
10008131	JNJ-10198409	55.63	15.325	12.125	3.89
10009557	SC-1	75.60	5.545	12.675	1.14

Figure 3. This heatmap shows the percentage of the total population within each of the indicated gates. Colors represent the deviation from the vehicle-treated cells (blue: within 3 SD of vehicle-treated cells; red: value of staurosporine-treated cells). Among the eleven flagged compounds, only staurosporine and PIK-75 (hydrochloride) were nearly completely cytotoxic, with most of the cell death being necrotic (>50% DAPI+). The remainder of the cytotoxic compounds induced highly mixed and variable populations, suggesting different mechanisms of cytotoxicity. For example, paclitaxel, JNJ-10198409, and 17β-hydroxy wortmannin were equally cytotoxic overall (around 60% live), but paclitaxel and JNJ-10198409 tended to induce more early and late apoptotic responses, while 17β-hydroxy wortmannin induced much more necrosis in these cells. In most cases, when TO-PRO[®]-3 was higher than average, annexin V was also higher than average, providing support that these independent markers measure related processes in apoptosis. However, in a few cases (e.g., NVP-BE2235), the percentage of cells which were TO-PRO[®]-3-positive was above average while the annexin V percentage was below average.

TMRE Staining of Treated Cells

Item No.	Item Name	Live	Early Apoptotic	Necrotic	Apoptotic	TMRE (% control)
10461	Paclitaxel	59.48	15.075	10.975	3.91	80.36
81590	Staurosporine	0.90	12.425	68.575	3.47	-0.62
10565	NVP-BE2235	74.73	8.295	11.1	1.08	154.19
11445	Tunicamycin	67.13	8.47	13.825	4.01	124.89
11491	AZD 7762	35.73	18.625	26.575	4.96	28.93
13812	17β-hydroxy Wortmannin	57.13	7.3525	29.975	0.69	125.94
11569	GSK1059615	76.18	6.4875	11.825	1.55	169.73
10009210	PIK-75 (hydrochloride)	2.29	20.425	55.775	10.83	0.86
10010248	AG-17	76.13	3.4775	16.4	0.47	0.62
10008131	JNJ-10198409	55.63	15.325	12.125	3.89	72.03
10009557	SC-1	75.60	5.545	12.675	1.14	108.32

Figure 4. The addition of TMRE staining (last column, percent of control-treated cells, blue; percentage of cells TMRE+ in vehicle, red; percentage of cells TMRE+ in staurosporine) to the apoptotic stains in this screen allows for the concomitant analysis of the health of the mitochondria, elucidating potential mitochondrial toxicity. While several compounds shown to induce apoptosis predictably had low levels of TMRE positivity, one stood out. Notably, AG-17 showed lower than average apoptosis induction, average necrosis levels, and almost no TMRE staining, suggesting a profound uncoupling of the mitochondria without direct induction of the apoptotic pathway. This conclusion is supported by literature describing mitochondrial disruption by AG-17 in the low micromolar range.¹

Conclusions

- Eleven compounds in the kinase screening library induced significant cell death.
- By analyzing the mode of cell death using the Early Apoptosis Detection Assay Kit (Cayman Item No. 601360), we were able to determine that there are differences in the way that each of these compounds kills, which may affect their therapeutic efficacy.
- We have shown, using the high-throughput MACSQuant[®] X Flow Cytometer by Miltenyi, that a library of compounds can be screened efficiently and accurately by flow cytometry.

Reference

1. Burger, A.M., Kaur, G., Alley, M.C., et al. Typhostin AG17, [(3S)-Di-tert-butyl-4-hydroxybenzylidene-malononitrile], inhibits cell growth by disrupting mitochondria. *Cancer Res.* 55(13), 2794-2799 (1995).

Learn what our Bioanalysis & Assay Development Services team can do for you
www.caymanchem.com/bioanalysis