




**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

# How to Choose a Microbiome Standard



## Controls and Standards in Microbiome Research

The advancement of NGS based technologies has led to a rapid growth in the field of microbiome research and deciphering microbial community composition, function, and interactions. Many studies conclude that technical variability in microbiome processing methods leads to significant variations in results<sup>1-3</sup>. Most of the discrepancies in reporting are explained by differences among the methods for nucleic acid extraction, NGS library preparation, bioinformatic data processing, and the choice of reference databases. Despite the complexity and variation introduced by varying protocols and methods for each step of the microbiomics workflow, data is being generated at an unprecedented pace. In many cases, a lack of proper controls or comparison to microbiome reference materials means that important and high-impact conclusions cannot be reproduced or reliably compared to similar data sets.

Commonly used and accepted controls or reference reagents are often called 'standards' because their inclusion and consideration allow for comparisons of methods, equipment, and protocols. Microbiome standards are imperative for microbial community profiling and analysis. Whereas the microbial compositions of experimental samples are variable and often unknown, microbiome standards provide a common, accurate, and consistent measurement as a basis for comparison. By providing a common control to measure and evaluate performance,

microbiome standards indicate biases allowing users to verify and optimize methods, enable inter-lab comparisons, and ensure reproducibility.

## How to Select the Appropriate Microbiome Controls

The principle of a microbiome standard is simple: use a well characterized, quantified, and known microbial input to perform experimental procedures and evaluate consistency of the output. Standards can then be run as a parallel quality control to experimental samples to evaluate the consistency of the method. The resulting profile provides a basis to calibrate and when needed, begin troubleshooting. Several different types of NGS Microbiome controls are available, each detecting different and sometimes overlapping parts of the complex microbiome processing workflow. This article is meant to aid in selecting the appropriate reference reagents and controls for your microbiome experiments.

## Mock Communities, True Diversity Reference, and Spike-in Controls

Several categories of microbiome reference reagents are available including microbial mock communities, true diversity reference material, and spike-in controls. Each category has overlapping characteristics, such as the use as positive controls, and each detects different biases throughout the microbiome analysis workflow.

Table 1 – Microbiome Standards and Controls Suggested Use

Mock Community Standards (Cellular)	
Standards	Suggested Applications
ZymoBIOMICS™ Microbial Community Standard	<ul style="list-style-type: none"> <li>General optimization and benchmarking</li> <li>Positive control for microbial lysis</li> </ul>
ZymoBIOMICS™ Gut Microbiome Standard	<ul style="list-style-type: none"> <li>General optimization and benchmarking for gut microbiome workflows</li> <li>Assess cross-kingdom, strain-level resolution, and pathogen detection</li> </ul>
ZymoBIOMICS™ Microbial Community Standard II (Log Distribution)	<ul style="list-style-type: none"> <li>Assessing detection limit of whole workflows beginning with DNA extractions</li> </ul>
Mock Community Standards (DNA)	
ZymoBIOMICS™ Microbial Community DNA Standard	<ul style="list-style-type: none"> <li>Optimization and positive control for library preparation and bioinformatics</li> </ul>
ZymoBIOMICS™ HMW DNA Standard	<ul style="list-style-type: none"> <li>Optimization and positive control for long-read sequencing library preparation and bioinformatics</li> </ul>
ZymoBIOMICS™ Microbial Community DNA Standard II (Log Distribution)	<ul style="list-style-type: none"> <li>Assessing detection limits of library preparation and bioinformatics</li> </ul>
True Diversity Reference	
ZymoBIOMICS™ Fecal Reference with TruMatrix™ Technology	<ul style="list-style-type: none"> <li>Assessing taxonomic assignment and bioinformatic processing parameters</li> <li>Enable inter-lab and inter-study data comparisons</li> </ul>
Spike-In Controls	
ZymoBIOMICS™ Spike-in Control I (High Microbial Load)	<ul style="list-style-type: none"> <li><i>In situ</i> extraction control and absolute quantification for high biomass samples</li> </ul>
ZymoBIOMICS™ Spike-in Control II (Low Microbial Load)	<ul style="list-style-type: none"> <li><i>In situ</i> extraction control and absolute quantification for low biomass samples</li> </ul>



distribution of species enables users to evaluate the detection limits of their microbiome analysis workflow<sup>8</sup>.

## DNA Mock Community Standards

Mock community standards made with purified microbial genomic DNA are more often used to detect biases and as optimization tools because they are utilized as input for library preparation rather than at the beginning of the workflow. DNA mock community standards such as the [ZymoBIOMICS™ Microbial Community DNA Standard](#) can be utilized to control biases associated with library prep and bioinformatics<sup>9-10</sup>. The optimization can be focused on library prep by first aligning NGS reads generated from the standard only to the genomes within the standard. After library prep has been optimized, the bioinformatics pipeline can be evaluated by aligning NGS reads against an entire reference database.

Similar to the cellular version, log distributed DNA standards, such as the [ZymoBIOMICS™ Microbial Community DNA Standard II \(Log Distribution\)](#), are used to assess detection limits but for library prep and bioinformatics pipelines.

Furthermore, an emerging technology for metagenomic analysis and genome assembly is long-read sequencing, often referred to as 3rd gen sequencing. Critical to long-read sequencing library prep and bioinformatics is high molecular weight DNA. The [ZymoBIOMICS™ HMW DNA Standard](#) is the only commercially available high molecular weight mock community, and has been used to evaluate sequencing chemistries and bioinformatic tools for long read sequencing<sup>11-12</sup>.

## True Diversity Reference

A true diversity reference is control material from a specified natural source that contains a complete, unchanging microbiome. In contrast to mock communities which have a quantified, known, and defined composition, the microbial composition of a true diversity reference is naturally derived. The [ZymoBIOMICS™ Fecal Reference with TruMatrix™ Technology\\*](#) is the first commercially available true diversity reference stabilized for long-term and lot-to-lot consistency. This reference features the high microbial diversity of a real fecal sample as well as a wide range of abundance.

Run-to-run and user-to-user consistency can be assessed on the same sample for each experiment. Reference materials can also be used to test system suitability by challenging experimental methods with actual source material. Bioinformatic analysis and taxonomy assignment are challenged with the added complexity of an unchanging true diversity sample. Since the microbial composition is static, the abundance and composition are stable and therefore allow users to assess method and analysis consistency.

## Spike-in Controls

Unlike mock communities and true diversity references, spike-in controls offer different functions when added directly to experimental samples. The [ZymoBIOMICS™ Spike-in Controls](#) are composed of very unique species, alien to the human microbiome as well as many others. This enables them to be spiked into samples without interfering with the native microbiome. The defined composition of these species enables the quantification of the absolute cell number within the unknown sample, when analyzed with NGS-based microbiome methods. Furthermore, an emerging use of these spike-in controls is as *in situ* quality controls, meaning that it can be used as a positive control for every sample rather than a positive control for a whole run. This is very useful for NGS-based pathogen diagnosis.

Two spike-in controls are available for different sample types. The [ZymoBIOMICS™ Spike-in Control I \(High Microbial Load\)](#) is meant for high biomass samples such as stool. The [ZymoBIOMICS™ Spike-in Control II \(Low Microbial Load\)](#) is meant for low microbial biomass samples such as sputum and bronchoalveolar lavage (BAL) fluid.

## Choosing a Microbiome Standard

The past several years has seen an explosion in the demand for microbiome standards, controls, and references that provide different and specific utilities. The scientists at Zymo Research share a passion for creating and providing the world with tools to improve microbiome data accuracy and reproducibility. As a result, the [ZymoBIOMICS™](#) line of standards, references, and controls provides a range of utility for various microbiome applications. Additional information about the standards and applications can be found in Table 2.

\*TruMatrix™ is a trademark of The BioCollective.

## References

1. Sinha R, Abu-Ali G, Vogtmann E, Fodor AA, Ren B, Amir A, Schwager E, Crabtree J, Ma S. Microbiome Quality Control Project C et al: Assessment of variation in microbial community amplicon sequencing by the Microbiome Quality Control (MBQC) project consortium. *Nat Biotechnol.* 2017; 35(11): 1077–86.
2. Costea PI, Zeller G, Sunagawa S, Pelletier E, Alberti A, Levenez F, Tramontano M, Driessen M, Hercog R, Jung FE, et al. Towards standards for human fecal sample processing in metagenomic studies. *Nat Biotechnol.* 2017; 35(11): 1069–76.
3. Jovel J, Patterson J, Wang W, Hotte N, O’Keefe S, Mitchel T, Perry T, Kao D, Mason AL, Madsen KL, et al. Characterization of the gut microbiome using 16S or shotgun metagenomics. *Frontiers in Microbiology.* 2016; 7:459.
4. Bartolomaeus TUP, Birkner T, Bartolomaeus H, Löber U, Avery EG, Mähler A, Weber D, Kochlik B, Balogh A, Wilck N, Boschmann M, Müller DN, Markó L, Forslund SK. Quantifying technical confounders in microbiome studies. *Cardiovascular Research.* 2021;17(3): 863-875.
5. Ojo-Okunola A, Claassen-Weitz S, Mwaikono KS, Gardner-Lubbe S, Zar HJ, Nicol MP, du Toit E. The Influence of DNA Extraction and Lipid Removal on Human Milk Bacterial Profiles. *MDPI Methods and Protocols.* 2020; 3(2): 39
6. Zhang B, Brock M, Arana C, Dende C, van Oers NS, Hooper LV, Raj P. Impact of bead-beating intensity on the genus and species level characterization of gut microbiome using amplicon and complete 16S rRNA gene sequencing. *Frontiers in Cellular and Infection Microbiology.* 2021; 11: 678522
7. Palkova L, Tomova A, Repiska G, Babinska K, Bokor B, Mikula I, Minarik G, Ostatnikova D, Soltys K. Evaluation of 16S rRNA primer sets for characterisation of microbiota in paediatric patients with autism spectrum disorder. *Nature Scientific Reports.* 2021; 11: 6781
8. Nicholls SM, Quick JC, Tang S, Loman NJ. Ultra-deep, long-read nanopore sequencing of mock microbial community standards. *GigaScience.* 2019; 8(5): giz043
9. Karst SM, Ziels RM, Kirkegaard RH, Sørensen EA, McDonald D, Zhu Q, Knight R, Albertsen M. High-accuracy long-read amplicon sequences using unique molecular identifiers with Nanopore or PacBio sequencing. *Nature Methods.* 2021; 18: 165-169.
10. Holm JB, Humphrys MS, Robinson CK, Settles ML, Ott S, Fu L, Yang H, Gajer P, He X, McComb E, Gravitt PE, Ghanem KG, Brotman RM, Ravel J. Ultrahigh-Throughput Multiplexing and Sequencing of >500-Base-Pair Amplicon Regions on the Illumina HiSeq 2500 Platform. *mSystems.* 2019; 4(1): e00029-19
11. Sereika M, Kirkegaard RH, Karst SM, Michaelsen TY, Sørensen EA, Wollenberg RD, Albertsen M. Oxford Nanopore R10.4 long-read sequencing enables near-perfect bacterial genomes from pure cultures and metagenomes without short-read or reference polishing. *bioRxiv.* 2021
12. Payne A, Holmes N, Clarke T, Munro R, Debebe BJ, Loose M. Readfish enables targeted nanopore sequencing of gigabase-sized genomes. *Nature Biotechnology.* 2021; 39: 442-450



**ZYMO RESEARCH**  
*The Beauty of Science is to Make Things Simple®*



[www.zymoresearch.com](http://www.zymoresearch.com)



[info@zymoresearch.com](mailto:info@zymoresearch.com)



Toll Free: (888) 882-9682

Learn More About the ZymoBIOMICS™ Microbiome Standards at

[www.zymoresearch.com/pages/microbiome-standards](http://www.zymoresearch.com/pages/microbiome-standards)