

---

Product Manual

# ViraDuctin™ AAV Transduction Kit


Catalog Number

AAV-200

10 transductions in 35 mm dish/6-well plate

**FOR RESEARCH USE ONLY**  
Not for use in diagnostic procedures

---



**CELL BIOLABS, INC.**

*Creating Solutions for Life Science Research*

## **Introduction**

The viral system includes vectors developed from retrovirus (RV), adenovirus (AdV), adeno-associated virus (AAV), lentivirus (LV), and herpes simplex virus (HSV). AAV belongs to the family of Parvoviridae, a group of viruses among the smallest of single-stranded and non-enveloped DNA viruses. There are eight different AAV serotypes reported to date.

Recombinant AAV-2 is the most common serotype used in gene delivery, and can be produced at high titers with a helper virus or Cell Biolabs' AAV Helper-Free System. AAV can infect both dividing and non-dividing cells and can be maintained in the human host cell, creating the potential for genome integration. Because AAV is a naturally defective virus, requiring provision of several factors in *trans* for productive infection, it is considered the safest viral vector to use. These factors make AAV an attractive vector for gene therapy.

The AAV transduction process includes viral binding and entry, intracellular trafficking, nuclear transport, and viral second strand DNA synthesis. The viral second strand DNA synthesis has been shown to be the rate limiting step, which leads to inefficient transduction by AAV vectors.

ViraDuctin™ AAV Transduction Kit is a proprietary formulation for the transduction of AAV. ViraDuctin™ provides the following advantages:

- Higher transduction efficiency in many cell types
- Easy to use
- Ideal for transduction of nonpermissive cells such as primary cells and stem cells

## **Related Products**

1. AAV-100: 293AAV Cell Line
2. VPK-140: ViraBind™ AAV Purification Kit
3. VPK-141: ViraBind™ AAV Purification Mega Kit
4. VPK-145: QuickTiter™ AAV Quantitation Kit
5. AD-200: ViraDuctin™ Adenovirus Transduction Reagent
6. LTV-200: ViraDuctin™ Lentivirus Transduction Kit
7. RV-200: ViraDuctin™ Retrovirus Transduction Kit

## **Kit Components**

1. ViraDuctin™ AAV Transduction Reagent A (20X) (Part No. 320007): One sterile tube – 1 mL.
2. ViraDuctin™ AAV Transduction Reagent B (100X) (Part No. 320008): One sterile tube – 200 µL.

## **Materials Not Supplied**

1. AAV Stock Solution
2. Cells and cell culture growth medium

### 3. 37°C Tissue Culture Incubator

#### **Storage**

Store kit components at -20°C. Avoid multiple freeze/thaws by aliquoting.

#### **Safety Considerations**

Remember that you will be working with samples containing infectious virus. Follow the recommended NIH guidelines for all materials containing BSL-2 organisms.

#### **Protocol**

The following transduction protocol is written for adherent cells in a 6-well plate or 35 mm culture dish. Refer to the below table for the appropriate dispensing volumes of other plate formats.

<b>Culture Dish</b>	<b>96-well</b>	<b>24-well</b>	<b>12-well</b>	<b>6-well or 35 mm</b>	<b>60-mm</b>	<b>10-cm</b>
ViraDuctin™ AAV Reagent A (20X) (μL)	5	25	50	100	250	500
ViraDuctin™ AAV Reagent B (100X) (μL)	1	5	10	20	50	100
Complete Culture Media (μL)	94	470	940	1880	4700	9400
Final Reagent Mixture (μL)	100	500	1000	2000	5000	10000



**Table 1: Dispensing Volumes of Different Plate Formats**

1. The day before transduction, trypsinize and count the cells, plating  $1-4 \times 10^5$  cells in 2.0-3.0 mL complete culture medium per well of a 6-well plate. Incubate cells at 37°C overnight to ensure firm attachment.
2. Warm the ViraDuctin™ AAV Transduction Reagents for at least 10 minutes at room temperature.
3. In a sterile tube, add 100 μL of ViraDuctin™ AAV Transduction Reagent A (20X) and 20 μL of ViraDuctin™ AAV Transduction Reagent B (100X), mix by inverting.
4. Incubate 5 minutes at room temperature.
5. Next, add 1.88 mL of complete culture medium to the ViraDuctin™ Reagent mixture, mix by inverting.
6. Incubate 5 minutes at room temperature.
7. Remove the culture medium from the overnight plated cells. Apply the entire volume of ViraDuctin™ Transduction Reagent mixture to the cells (~ 2 mL).
8. Incubate at 37°C overnight.
9. Carefully, aspirate the media from the cells.
10. Wash the cells twice with 2 mL of complete culture medium.
11. Add 1 mL of desired AAV-containing media to the plate well. Incubate at 37°C for 1-2 hours, gently swirling/mixing every 30 minutes.

12. Add an additional 1 mL of pre-warmed complete media per well, incubating at 37°C for 2-7 days.

### Example of Results

The following figures demonstrate typical transduction results. One should use the data below for reference only. This data should not be used to interpret actual results.



**Figure 1: AAV-GFP Transduction Efficiency with ViraDuctin™.** HeLa cells were seeded at 50,000 cells/well in a 24-well plate overnight. Cells were then treated with ViraDuctin™ for 20 hours before infecting with AAV-GFP.

## **References**

1. Rabinowitz, J, and Samulski, R. J. (1998). *Curr. Opin. Biotechnol.*, **9**, 470-475.
2. Summer ford, C., and Samulski, R. J. (1999). *Nat. Med.*, **5**, 587-588.
3. Clark, K., Liu, X., McGrath, J., and Johnson, P. (1999). *Hum. Gene Ther.*, **10**, 1031-1039.

## **Recent Product Citations**

1. Lau, C-H. et al. (2019). Targeted Transgene Activation in the Brain Tissue by Systemic Delivery of Engineered AAV1 Expressing CRISPRa. *Molecular Therapy: Nucleic Acid*. doi: 10.1016/j.omtn.2019.04.015.
2. Nazari, M. et al. (2014). AAV2-mediated follistatin overexpression induces ovine primary myoblasts proliferation. *BMC Biotechnol.* **14**:87.

## **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

## **Contact Information**

Cell Biolabs, Inc.  
7758 Arjons Drive  
San Diego, CA 92126  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2008-2019: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.