
Product Manual


CytoSelect™ 96-Well Phagocytosis Assay (Zymosan Substrate)

Catalog Number

CBA-224

96 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research

Introduction

In mammals, phagocytosis by phagocytes (e.g., macrophages, dendritic cells, and neutrophils) is essential for a variety of biological events, including tissue remodeling and the continuous clearance of dying cells. Furthermore, phagocytosis represents an early and crucial event in triggering host defenses against invading pathogens. Phagocytosis comprises a series of events, starting with the binding and recognition of particles by cell surface receptors, followed by the formation of actin-rich membrane extensions around the particle. Fusion of the membrane extensions results in phagosome formation, which precedes phagosome maturation into a phagolysosome. Pathogens inside the phagolysosome are destroyed by lowered pH, hydrolysis, and radical attack (Figure 1). These early events that are mediated by the innate immune system are critical for host survival. As a result of this process, pathogen-derived molecules can be presented at the cell surface (antigen presentation), allowing the induction of acquired immunity.





Figure 1: Phagocytosis Processes.

Zymosan (*Saccharomyces cerevisiae*) is prepared from yeast cell wall and consists of protein-carbohydrate complexes. Zymosan is a commonly used pathogen in phagocytosis assays. Typically, engulfed Zymosan particles are manually counted (expressed as a phagocytosis index or engulfed particles per phagocyte). This manual counting method is quite cumbersome, time-consuming, and difficult when testing a large number of samples.

Cell Biolabs' CytoSelect™ 96-well Phagocytosis Assay (Zymosan) uses prelabeled Zymosan particles as a phagocytosis pathogen; however, it does **not** involve subjective manual counting of Zymosan particles inside cells. Instead external Zymosan particles are blocked before the colorimetric detection of engulfed particles (Figure 2). This format provides a quantitative, high-throughput method to

accurately measure phagocytosis. The CytoSelect™ 96-well Phagocytosis Assay (Zymosan) provides a robust system for screening TLR ligands, phagocytosis activators or inhibitors. Each kit provides sufficient quantities to perform 96, 48, 24 tests in a 96, 48, 24-well plate, respectively.

Assay Principle



Related Products

1. CBA-220: CytoSelect™ 96-Well Phagocytosis Assay (Red Blood Cell, Colorimetric Format)
2. CBA-210: CytoSelect™ Leukocyte-Endothelium Adhesion Assay
3. CBA-211: CytoSelect™ Leukocyte-Epithelium Adhesion Assay
4. CBA-212: CytoSelect™ Leukocyte Transmigration Assay

Kit Components

1. Zymosan Suspension (Part No. 122401): One 1 mL tube of prelabeled Zymosan in PBS, 5×10^8 particles/mL
2. Fixation Solution (Part No. 122402): One 20 mL bottle of 3.2% Buffered Formaldehyde Solution
3. Blocking Reagent (100X) (Part No. 122403): One 200 μ L tube
4. 10X Permeabilization Solution ((Part No. 122404): One 1.5 mL tube of PBS/1% Triton X-100
5. Detection Reagent (250X) (Part No. 122405): One 50 μ L tube
6. Detection Buffer (Part No. 122406): One 10 mL bottle
7. Substrate (Part No. 122407): One 12 mL amber bottle
8. Stop Solution (Part No. 122408): One 12 mL bottle
9. Phagocytosis Inhibitor (Part No. 122006): One amber tube – 20 μ L of 2 mM Cytochalasin D in DMSO

Materials Not Supplied

1. Phagocyte and Culture Medium
2. PBS, PBS/0.1% BSA
3. 37°C Incubator, 5% CO₂ Atmosphere
4. Light Microscope
5. 96-well Mitrotiter Plate
6. 96-well Microtiter Plate Reader

Storage

Store all kit components at 4°C.

Preparation of Reagents

- **Zymosan Suspension:** Thaw Zymosan suspension at 4°C. Either nonopsonized or opsonized Zymosan particles can be used in phagocytosis assay. To opsonize Zymosan particles, incubate particles with desired serum or IgG for 30 minutes at 37°C, pellet particles by centrifugation and wash a few times with sterile 1X PBS. Prior to using, resuspend the particles in the same volume of sterile 1X PBS. Store at 4°C.
- **1X Blocking Reagent:** Prepare the appropriate volume for the number of samples being tested. IMMEDIATELY prior to using, dilute the provided 100X Blocking Reagent 1:100 in 1X PBS/0.1% BSA. Do not store.
- **1X Permeabilization Solution:** Prepare the appropriate volume for the number of samples being tested. Prior to using, dilute the provided 10X Permeabilization Solution 1:10 in 1X PBS. Store at 4°C.

- 1X Detection Reagent: Prepare the appropriate volume for the number of samples being tested. IMMEDIATELY prior to using, dilute the provided 250X Detection Reagent 1:250 in 1X PBS/0.1% BSA. Do not store.

Assay Protocol: Adherent Phagocytes

The following assay protocol is written for a 96-well format. Refer to the below table for the appropriate dispensing volumes of other plate formats.

Culture Dish	96-well	48-well	24-well
Phagocyte Seeding Volume (μL/well)	100	200	400
Zymosan Suspension (μL/well)	10	20	40
Fix Solution (μL/well)	100	200	400
Permeabilization Solution (μL/well)	100	200	400
Detection Buffer (μL/well)	50	100	200

Table 1: Dispensing Volumes of Different Plate Formats

I. Phagocytosis of Zymosan

1. Harvest and resuspend phagocytic cells in culture medium at $1 - 5 \times 10^5$ cells/mL or the appropriate concentration that yields 50-80% confluency after overnight incubation. Seed 100 μL in each well of a 96-well plate and incubate overnight at 37°C, 5% CO₂.
2. Treat phagocytes with desired activators or inhibitors.
3. Add 10 μL of Zymosan suspension to each well. Mix well and immediately transfer the plate to a cell culture incubator for 15 minutes – 2 hours. Each sample including a negative control without Zymosan particles should be assayed in duplicate.
4. Remove the culture medium by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times.
5. Gently add 200 μL of cold, serum-free medium (e.g. DMEM, RPMI) to each well. Promptly remove the cold media by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times. Repeat twice more.

Note: For loosely attached cells, complete culture media is preferred to maintain cell attachment

II. Remove and block external particles

Note: Perform steps with care, gently adding solutions as to not disrupt cell attachment

1. Add 100 μL of Fixation Solution to each well, incubating 5 minutes at room temperature.
2. Promptly remove the Fixation Solution by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times.

3. Wash twice with 1X PBS.
4. Add 100 μ L of prediluted 1X Blocking Reagent to each well (see Preparation of Reagents Section). Incubate the plate for 60 minutes at room temperature on an orbital shaker.
5. Promptly remove the Blocking Reagent by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times. Wash three times with 1X PBS.

III. Detection of internalized particles

Note: Perform steps with care, gently adding solutions as to not disrupt cell attachment

1. Remove the PBS wash and add 100 μ L of prediluted 1X Permeabilization Solution (see Preparation of Reagents Section) to each well, incubate 5 minutes at room temperature.
2. Promptly remove the Permeabilization Solution by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times. Wash once with 1X PBS.
3. Add 100 μ L of prediluted 1X Detection Reagent to each well (see Preparation of Reagents Section). Incubate the plate for 60 minutes at room temperature on an orbital shaker.
4. Promptly remove the Detection Reagent Solution by gently aspirating or inverting the plate and blotting on a paper towel. Gently tap several times. Wash three times with 1X PBS.
5. Add 50 μ L of Detection Buffer to each well. Incubate the plate for 10 minutes at room temperature on an orbital shaker.
6. Initiate the reaction by adding 100 μ L of Substrate. Incubate for 5-20 minutes at 37°C.
7. Stop the reaction by adding 50 μ L of the Stop Solution and mix by placing the plate on an orbital plate shaker for 30 seconds.
8. Read the absorbance of each well at 405 nm.

Assay Protocol: Suspension Phagocytes

1. Harvest and resuspend phagocytic cells in culture medium at $0.2 - 1.0 \times 10^6$ cells/mL. Seed 100 μ L in each well of a 96-well plate.
2. Treat phagocytes with desired activators or inhibitors.
3. Add 10 μ L of Zymosan suspension to each well. Mix well and immediately transfer the plate to a cell culture incubator for 15 minutes – 2 hours.
4. Remove the culture medium by centrifugation and gentle aspiration.
5. Add 200 μ L of cold 1X PBS to each well. Promptly remove the PBS Solution by centrifugation and gentle aspiration.
6. Add 100 μ L of Fixation Solution to each well, incubate 5 minutes at room temperature.
7. Promptly remove the Fixation Solution by centrifugation and gentle aspiration.
8. Wash twice with 1X PBS.
9. Add 100 μ L of prediluted 1X Blocking Solution to each well (see Preparation of Reagents Section). Incubate the plate for 60 minutes at room temperature on an orbital shaker.

10. Promptly remove the Blocking Solution by centrifugation and gentle aspiration. Wash three times with 1X PBS.
11. Add 100 μL of prediluted 1X Permeabilization Solution (see Preparation of Reagents Section) to each well, incubate 5 minutes at room temperature.
12. Promptly remove the Permeabilization Solution by centrifugation and gentle aspiration. Wash twice with 1X PBS.
13. Add 100 μL of prediluted 1X Detection Reagent to each well (see Preparation of Reagents Section). Incubate the plate for 60 minutes at room temperature on an orbital shaker.
14. Promptly remove the Detection Reagent Solution by centrifugation and gentle aspiration. Wash three times with 1X PBS.
15. Add 50 μL of Detection Buffer to each well and incubate the plate for 10 minutes at room temperature on an orbital shaker.
16. Initiate the reaction by adding 100 μL of Substrate. Incubate for 5-20 minutes at 37°C.
17. Stop the reaction by adding 50 μL of the Stop Solution and mix by placing the plate on an orbital plate shaker for 30 seconds.
18. Read the absorbance of each well at 405 nm.

Example of Results

The following figures demonstrate typical results with the CytoSelect™ 96-well Phagocytosis Assay Kit. Absorbance measurements were performed on a Microplate Autoreader EL311 (Bio-Tek Instruments Inc.) with a 405 nm filter. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 2. Inhibition of Raw 264.7 Macrophage Phagocytosis by Cytochalasin D. 50,000 cells/well of Raw 264.7 macrophages were seeded overnight in a 96-well plate. Cytochalasin D was used to pretreat Raw 264.7 cells for 1 hr at 37°C before addition of Zymosan particles at 50:1 ratio. Phagocytosis was stopped after 30 minutes and the amount of engulfed Zymosan particles was determined as described in the Assay Protocol.




Figure 3. Zymosan Particles Engulfment by Raw 264.7 Macrophage.

References

1. Sansonetti, P. (2001) *Semin. Immunol.* **13**:381–390.
2. Jutras I and Desjardins M. (2005) *Annu Rev Cell Dev Biol.* **21**:511-27.
3. Janeway, C. A., Jr., and Medzhitov R. (2002) *Annu. Rev. Immunol.* **20**:197–216.
4. Gordon, S. (2002) *Cell.* **111**:927–930.

Recent Product Citations

1. Lee, S. G. et al. (2015). Immunostimulatory polysaccharide isolated from the leaves of *Diospyros kaki* Thumb modulate macrophage via TLR2. *Int J Biol Macromol.* **79**:971-982.
2. Jung, J. Y. et al. (2015). *Lactobacillus sakei* K040706 evokes immunostimulatory effects on macrophages through TLR 2-mediated activation. *Int Immunopharmacol.* doi:10.1016/j.intimp.2015.05.037.
3. Zhang, H. et al. (2015). Functional analysis and transcriptomic profiling of iPSC-derived macrophages and their application in modeling mendelian disease. *Circ Res.* doi:10.1161/CIRCRESAHA.117.305860.
4. Fiorcari, S. et al. (2015). Lenalidomide interferes with tumor-promoting properties of nurse-like cells in chronic lymphocytic leukemia. *Haematologica.* **100**:253-262.
5. Liao, W. T. et al. (2014). Cyclic GMP-dependent protein kinase II is necessary for macrophage M1 polarization and phagocytosis via toll-like receptor 2. *J Mol Med (Berl).* doi: 10.1007/s00109-014-1236-0.
6. Kasat, K. et al. (2014). Anti-inflammatory actions of endogenous and exogenous interleukin-10 versus glucocorticoids on macrophage functions of the newly born. *J Perinatol.* **34**:380-385.

7. Haselow, K. et al. (2013). Bile acids PKA-dependently induce a switch of the IL-10/IL-12 ratio and reduce proinflammatory capability of human macrophages. *J. Leukoc. Biol.* **94**:1253-1264.
8. Pierce, L.M. et al. (2012). Effect of heavy metal tungsten alloy particles on oxidative product formation and phagocytosis by lung macrophages. *Am. J. Respir. Crit. Care Med.* **185**:A4666.
9. Polancec, D.S. et al. (2012). Azithromycin drives in vitro GM-CSF/IL-4-induced differentiation of human blood monocytes toward dendritic-like cells with regulatory properties. *J Leukoc Biol.* **91**:229-243.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2007-2015: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.