


CALCIUM INDICATORS SELECTION GUIDE

Cat no.	Indicator Dye	Ex	Em	Laser	Filter	Kd
20412	Cal-500™	388	482	405	DAPI	303 nM
20650	Calbryte™ 520	493	515	488	FITC	1.2 μM
21130	Cal-520®	493	515	488	FITC	320 nM
20642	Calbryte™-520L	493	515	488	FITC	91 μM
21142	Cal-520FF™	493	515	488	FITC	9.8 μM
21146	Cal-520N™	493	515	488	FITC	90 μM
20507	OG488 BAPTA -1	493	522	488	FITC	170 nM
21080	Fluo-8®	495	516	488	FITC	389 nM
21090	Fluo-8H™	495	516	488	FITC	232 nM
21096	Fluo-8L™	495	516	488	FITC	1.9 μM
21104	Fluo-8FF™	495	516	488	FITC	10 μM
20550	Fluo-4	495	528	488	FITC	345 nM
20501	Cal Green™ 1	498	517	488	FITC	190 nM
21011	Fluo-3	506	515	488	FITC	390 nM
21120	Rhod-4™	523	551	561	TRITC	451 nM
21062	Rhod-2	553	577	561	TRITC	570 nM
21077	Rhod-FF	553	577	561	TRITC	19 μM
21070	Rhod-5N	557	580	561	TRITC	300 μM
20510	Cal-590™	574	588	561	Cy3/TRITC	561 nM
20700	Calbryte™ 590	581	593	561	Cy3/TRITC	1.4 μM
20720	Calbryte™ 630	607	624	561	Texas Red	1.2 μM
20530	Cal-630™	609	626	561	Texas Red	792 nM
20455	Cal-670™	667	680	640	Cy5	853 nM
20460	Cal-770™	758	783	750	Cy7	850 nM

* FCa/FFree = Increase in fluorescence intensity of calcium bound indicator relative to its calcium free indicator form.



◀ Response of endogenous P2Y receptor to ATP in CHO-K cells. CHO-K cells were seeded overnight at 40,000 cells per 100 μL per well in a 96-well black wall/clear bottom costar plate. 100 μl of 4 μM Cal 520™ AM (cat. no. 21130) in HHBS with 1 mM probenecid were added into the wells, and the cells were incubated at 37 °C for 1 hour. The dye loading mediums were replaced with 100 μl HHBS and 1 mM probenecid, then imaged with a fluorescence microscope (Olympus IX71) using FITC channel before and after adding 50 μl of 300 μM ATP.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. © 2025 AAT Bioquest, Inc. All rights reserved. All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of AAT Bioquest.

