

THE SWIFT RNA LIBRARY KIT OPTIMIZES RNA-SEQ DATA QUALITY AND COSTS FOR FFPE SAMPLES

SUMMARY

This application note illustrates how the Swift RNA Library Kit enables multiple workflow options to maximize the output of RNA-Seq data from FFPE samples, including upstream ribodepletion of total RNA and downstream enrichment of total RNA libraries by hybridization capture. For the analysis of coding sequence information, a 4-fold savings in sequencing cost per sample can be achieved using the hybridization capture workflow. The information and methods presented below allow core facilities, production labs, and research scientists to choose the best workflow for their sequencing goals and to reduce costs and time associated with RNA-Seq library preparation and sequencing.

- For FFPE samples, DV₂₀₀ is more informative than RIN score
- Ribodepletion preserves coding and non-coding RNA but is not cost-effective for coding sequence analysis of FFPE samples due to only 20-30% of reads representing exonic regions
- Hybridization capture provides the most cost-effective workflow for coding sequence analysis, with 6-plex library capture resulting in a 2-fold savings for FFPE RNA-to-fastq costs
- Comparative analysis of both workflows details the RNA recovered, kits and reagents used, preparation time, and costs per sample

INTRODUCTION

The majority of today's Next Generation Sequencing (NGS) libraries are made from RNA samples for the quantitative evaluation of biological function. Among the most prevalent RNA-Seq samples are RNA isolated from human tissues preserved through formalin fixation and paraffin-embedding (FFPE). Although this practice is effective for long-term sample preservation, it introduces several challenges for NGS, such as extracting high-integrity RNA, achieving sufficient library yield, and generating quality data. Commonly utilized commercial and homebrew techniques for RNA-Seq library preparation from FFPE samples are vulnerable to crosslinked or fragmented RNA substrates as well as to an increased representation of intronic reads in the final library, an observation that is consistently seen for formalin-fixed samples¹. In addition, libraries produced from limited material suffer from high duplication rates and adapter dimers, as well as lower mapping rates and fewer genes and transcripts detected.

As a solution for FFPE samples, we present the Swift RNA Library Kit. This stranded RNA-Seq kit leverages patented Adaptase[®] technology, a method for adapter ligation to single-stranded substrates, to directly convert first-strand cDNA into NGS libraries in 4.5 hours (**Figure 1**). This approach saves time, guarantees stranded libraries without requiring conventional second-strand cDNA synthesis or degradation, and maximizes yield and complexity while minimizing bias and undesired library

byproducts. Further, Swift RNA can make quality libraries from as little as 10 ng total FFPE RNA, allowing it to accommodate samples that have a limited starting input.

Most labs consider the depletion of ribosomal RNA (rRNA) (known as ribodepletion) as the default approach for processing FFPE RNA samples – despite high costs for reagents and sequencing – since poly(A)-selection would only provide 3' coverage of fragmented transcripts.




Figure 1: Swift RNA Library Kit Workflow

However, hybridization capture can also be used to produce quality data from FFPE RNA-Seq libraries. The Swift RNA kit is compatible with both of these FFPE workflows (**Figure 2**). The choice between ribodepletion and hybridization capture should be informed by the RNA types and sequencing data of interest. Below, we discuss each workflow and provide representative data from

Universal Human Reference (UHR) RNA and multiple breast cancer tumor FFPE samples. Additionally, we provide a comparison of sample types, workflow time, and library and sequencing costs to help inform decisions made for maximizing recovery and data quality while minimizing costs.




Figure 2: RNA-Seq Library Workflows - Upstream Ribodepletion vs. Downstream Enrichment

RNA CHARACTERIZATION

The quality of input RNA into an RNA-Seq workflow can substantially impact the library yield and resulting data quality. Because FFPE RNA is often degraded, assessing the starting quality is important for adjusting the RNA-Seq workflow to optimize results. RNA quality and integrity can be evaluated using two different metrics: RIN score (RNA Integrity Number; ratio of the 28S to 18S rRNA peaks) and DV₂₀₀ (percentage of RNA fragments that are greater than 200 nucleotides).

RIN scores are useful for determining whether a poly(A)-selection module can be used, as they are typically only recommended for samples with RIN ≥ 7 . High-quality samples typically have RIN scores ≥ 7 whereas FFPE samples often have RIN scores < 7 (**Figure 3**). For samples with RIN < 7 , ribodepletion and hybridization capture are the preferred workflows. For these samples, more information about sample quality can be obtained through the DV₂₀₀ score which can help inform the fragmentation time, SPRI ratio, and PCR cycling. In general, it is not recommended to make libraries using samples with a RIN < 2 or a DV₂₀₀ < 30 .

We obtained FFPE curls from breast cancer tumors (Spectrum Health; Grand Rapids, MI) and extracted the RNA using the RNeasy FFPE Kit (Qiagen 73504). Isolated RNA was analyzed using an RNA 6000 Pico Kit (Agilent 5067-1513) on the Agilent Bioanalyzer. Trace

analyses show all samples have a RIN score near 2 whereas the DV₂₀₀ scores range from 48 to 73 (**Figure 3**), providing a more accurate depiction of RNA integrity.




Figure 3: RNA Bioanalyzer Traces for UHR and FFPE RNA

RIBODEPLETION WORKFLOW

Ribodepletion of FFPE RNA is the recommended workflow when coding and noncoding sequence information is desired. For intact RNA samples, poly(A) selection is considered a simpler and lower-cost approach to enrich the sample for mRNA. However, because FFPE samples are damaged, mRNA transcripts may be broken away from their poly(A) tails, preventing this method from enriching for transcripts in their entirety. Ribodepletion, however, can work effectively on FFPE samples even if they are damaged because it depletes rRNA rather than selects for mRNA, therefore enabling the recovery of fragments without poly(A) tails. Additionally, because ribodepletion leaves behind all non-rRNA molecules, it is the preferred workflow for those interested in intronic regions, such as those from nuclear pre-mRNA and long non-coding RNAs.

Here, we present data from UHR RNA and two FFPE RNA samples with variable DV₂₀₀ scores (Samples 1 and 2, **Figure 3**) that underwent ribodepletion prior to library prep using the Swift RNA Library Kit (Cat. No. R1024/96) in duplicate. Using 100 ng of total RNA, all samples showed robust library yields with high mapping rates and preservation of transcript expression, as evidenced by a high correlation of the ERCC spike-in² transcripts (**Table 1**). Further, FFPE libraries detected an ample and

Sample	Library yield (nM)	STAR mapping rate (%)	Genes detected	Exonic rate (%)	rRNA rate (%)	ERCC (R ²)
UHR RNA	44.2	86.2	17,132	46.0	3.4	0.94
	33.6	87.5	17,180	46.2	3.2	0.92
FFPE Sample 1	14.0	83.6	15,095	22.0	6.7	0.91
	19.8	85.7	15,112	21.2	6.0	0.92
FFPE Sample 2	24.9	84.2	15,274	24.6	4.8	0.92
	23.9	85.5	15,325	25.7	4.5	0.94

Table 1: RNA-Seq Data for Ribodepleted FFPE

Ribodepletion Workflow Methods

50-100 ng of UHR and FFPE Samples 1-2 were ribodepleted using Lexogen Ribocop⁴ V1.3 and eluted in 5 μ L dH₂O for input into the Swift RNA Library Kit⁵. Swift RNA libraries were made with the following adjustments: fragmentation at 65°C for 5 min and a 1.2X post-RT SPRI clean-up (see **Protocol Considerations** below). All libraries were amplified using 12 PCR cycles. Libraries were sequenced on a Miniseq 2 x 75 bp and downsampled to 3 million reads prior to analysis with STAR⁶, picard⁷, rnaseq⁸, and fgbio⁹. CDS = coding sequence; UTR = untranslated region (5' and 3').

HYBRIDIZATION CAPTURE WORKFLOW

Hybridization capture of FFPE RNA-Seq libraries enables the capture and sequencing of specific regions of the transcriptome. For example, Whole Exome panels capture only exonic fragments, maximizing the representation of coding regions in the sequencing data. Alternately, projects with narrower focus can use panels designed to target cancer (Swift Cat. No. 83316) or inherited disease (Swift Cat. No. 83416) genes, among others. Although hybridization capture is a longer workflow compared to ribodepletion, it can dramatically reduce overall costs by capturing and sequencing only regions of interest (see **Workflow Tradeoffs** below).

Library Prep Adjustments for Ribodepleted FFPE RNA:

- Fragmentation time
- Post-RT SPRI Clean-up

See "Protocol Considerations" on Page 5.

consistent number of genes, despite variability in DV₂₀₀ scores (**Table 1**). More genes were detected for UHR RNA, likely due to increased RNA integrity, higher sample complexity, and a higher proportion of exonic reads.

Exonic rates are typically low for ribodepleted RNA (**Table 1**). Because ribodepletion only removes rRNA, it allows intronic and intergenic RNA to become library molecules, thereby reducing the representation of exonic fragments (40-50% for UHR RNA). The exonic rate for FFPE samples is even lower due to the bias toward intronic reads that comes as a byproduct of formalin-fixation¹ as well as due to tissue-specific profiles, such as an increased prevalence of intronic reads in breast cancer tumors³ (**Figure 4**).




Figure 4: Read Distribution for Ribodepleted RNA-Seq Libraries

Library Prep Adjustments for Total FFPE RNA meant for Hybridization Capture:

- Fragmentation time
- All SPRI Clean-ups
- PCR cycles

See "Protocol Considerations" on Page 5.

Here, we present data from UHR and FFPE RNA (Samples 1 and 2, **Figure 3**) processed using the Swift RNA Library Kit with hybridization capture adjustments (see **Protocol Considerations** below). Libraries were enriched for human coding sequences using the Swift Hyb, Wash, and Blocker kit (Cat. No. 89016) and Swift Exome Hyb Panel¹⁰

(Cat. No. 83216). The Exome Panel contains over 400,000 probes and spans a 39 Mb target region comprising 19,396 genes. The use of probes designed specifically to the exome more than triples the exonic rate compared to ribodepleted samples (>80% compared to ~25%, respectively; see **Tables 1-2** and **Figure 5**). Additionally, hybridization capture provides consistent data metrics for both intact (UHR) and damaged (FFPE) samples with a similar number of genes being detected (**Table 2**). Note that **Table 2** data is at a sequencing depth of 3 million reads per sample; deeper sequencing captures all 19,396 genes targeted by the Exome Panel.

Sample	Library yield (ng/uL)	STAR mapping rate (%)	Genes detected	Exonic rate (%)	rRNA rate (%)	Selected bases (%)
UHR RNA	82	89.6	17,758	91.2	1.3	96.7
FFPE Sample 1	93	88.8	17,009	89.2	0.9	96.5
FFPE Sample 2	73	87.9	17,024	88.1	1.0	95.8

Table 2: RNA-Seq Data for Hybridization Capture (100 ng input)

Robust Data Across FFPE Samples

The variability in sample quality, particularly for FFPE samples, often acts as a hurdle to producing consistent, quality libraries. In order to evaluate a wider range of samples, we performed hybridization capture on Swift RNA libraries from four additional breast cancer tumor FFPE samples with DV₂₀₀ scores ranging from 54 to 62 (**Figure 3**). Because FFPE samples are often limited in yield following RNA extraction, we also tested the ability of the workflow to support an input of 10 ng FFPE RNA.

Swift RNA Libraries followed by hybridization capture (Swift Exome Panel) result in sequencing data with high yields, mapping rates, and exonic rates across a range of FFPE samples (**Table 3**). Further, libraries also produced high-quality data despite the lower input of 10 ng.

Sample	Library yield (ng/uL)	STAR mapping rate (%)	Genes detected	Exonic rate (%)	rRNA rate (%)	Selected bases (%)
FFPE Sample 3	23.0	89.4	16,089	87.4	1.5	95.1
FFPE Sample 4	54.0	88.0	16,065	88.7	1.3	95.6
FFPE Sample 5	33.4	87.5	15,653	87.6	1.3	94.6
FFPE Sample 6	56.0	87.6	16,246	88.4	1.2	95.4

Table 3: RNA-Seq Data for Hybridization Capture (10 ng input)

Hybridization Capture Methods

Either 100 ng (**Table 2**) or 10 ng (**Table 3**) of the respective FFPE Samples (DV₂₀₀ scores from 48 to 73, see **Figure 3**) were used as input into the Swift RNA Library Kit⁵. Libraries were made with fragmentation at 94°C for 15 min (UHR) or 2 min (FFPE). All SPRI clean-ups in the protocol were adjusted to 1.8X. 100 ng libraries used 11 pre-hyb PCR cycles and 10 ng libraries used 15 pre-hyb PCR cycles before input into the Swift Exome Hyb Panel¹⁰. Pooled samples used either 9 (2-plex, **Figure 6**), 8 (4-plex, **Tables 2-3**), or 7 (2-plex, **Figure 6**) post-hyb PCR cycles. Libraries were sequenced on a Miniseq 2 x 75 bp and downsampled to 3 million reads prior to analysis with STAR⁶, picard⁷, and rnaseq⁸. CDS = coding sequence; UTR = untranslated region (5' and 3').




Figure 5: Read Distribution for Hybridization Capture of RNA-Seq Libraries

Multiplex up to 6 Libraries

RNA libraries can be combined prior to hybridization capture to maximize the value of each reaction. For DNA libraries, samples can often be multiplexed with up to 16 samples. However, because transcripts are expressed at variable levels, RNA libraries run the risk of depleting probes for highly-expressed genes, thus limiting the number of samples that can be multiplexed per capture.

Here, we show FPKM (Fragments Per Kilobase of transcript per Million mapped reads) plots for UHR and FFPE RNA libraries that have been multiplexed with 2, 4, or 6 samples. Linear regression of the same library with different multiplexing shows extremely high correlation (**Figure 6**). Multiplexing can reduce the overall cost of each sample (See **Workflow Tradeoffs** below).




Figure 6: FPKM Correlation of Multiplexed Hyb Capture Libraries

WORKFLOW TRADEOFFS

Two workflows were presented above for processing FFPE samples: 1) ribodepletion, which involves the removal of rRNA prior to library prep and 2) hybridization capture, which involves enriching exonic library molecules (e.g. using a Whole Exome panel) from total RNA libraries. The table below summarizes these two workflows, including the types of RNA recovered, reagents required, workflow details, and costs.

Table 4: Workflow and Cost Comparison for Ribodepletion vs. Hybridization Capture

WORKFLOW	RIBODEPLETION-BASED	HYB-CAPTURE BASED
RNA Types Recovered	Coding and non-coding RNA	Coding RNA
Kits Used	<ol style="list-style-type: none"> 1. Lexogen RiboCop V1.3 H/M/R (037.96) 2. Swift RNA Library Kit (R1096) 3. Swift Unique Dual Indexing Primer Kit (X9096) 	<ol style="list-style-type: none"> 1. Swift RNA Library Kit (R1096) 2. Swift Unique Dual Indexing Primer Kit (X9096) 3. Swift Exome Hyb Panel (83216) 4. Swift Hyb, Wash, and Universal Blocker Kit (89016) 5. Swift Library Amplification Primer Mix (88196) and KAPA 2X HiFi HotStart ReadyMix (KK2602)
Number of Steps	20	28
Library Prep Time (RNA-to-Library)	6 hours	12 hours
Library Price Per Sample (RNA-to-Library)	\$61	\$104 (6-plex libraries per capture)
Sequencing Price Per Sample Novaseq SP 100 cycle (Library-to-Fastq, 30M exonic reads)	\$472	\$121
Total Price Per Sample (FFPE RNA-to-Data, \$100/hr FTE)	\$566	\$284

**Estimates are based on US 2020 List prices and assume 24 samples (4 captures) are being prepared simultaneously.*

PROTOCOL CONSIDERATIONS

FFPE Samples are often degraded and require workflow-specific adjustments to the protocol in order to optimize the yield and quality of the final libraries. Follow the guidelines below and see [Table 5](#) for more information.

Fragmentation Time

RNA is fragmented following a chemical- and heat-based process. The amount of fragmentation can be adjusted by increasing or decreasing the fragmentation time and/or temperature. Because FFPE RNA is often degraded, it is already fragmented to a certain extent. The DV₂₀₀ score can help determine how fragmented the RNA is. For samples with low DV₂₀₀ scores, reduce the temperature and time when following the ribodepletion workflow. When following the hybridization capture workflow, reduce only the time (see [Table 5](#)). These recommendations can be further adjusted to optimize the fragmentation based on specific samples or insert size requirements.

SPRI Clean-up Ratios

SPRI clean-ups act as a size selection mechanism as well as a buffer exchange. A more relaxed SPRI (i.e., a

higher ratio of beads to sample volume) will result in carry over of smaller fragments whereas a more stringent SPRI (i.e., a lower ratio of beads to sample volume) will select for larger fragment sizes. Because FFPE RNA is often degraded, relaxing the SPRI ratio after reverse transcription (ribodepletion workflow) or throughout the protocol (hybridization capture workflow) preserves smaller fragments that can then be converted into library molecules (see [Table 5](#)).

A more relaxed SPRI ratio, however, can also enable the carry-over of adapter dimers, library artefacts that arise when two adapter molecules become ligated together. Adapter-dimers can cluster during sequencing to take up valuable sequencing reads and are often exacerbated at lower inputs. For the ribodepletion workflow, if adapter dimers are present in the final libraries, as evidenced by a peak near ~150 bp on an electrophoretic instrument, a second post-PCR SPRI can be performed. For the hybridization capture workflow, adapter dimers do not present an issue as they will not be captured during the hybridization steps. Thus, libraries that are meant for hybridization capture can utilize relaxed 1.8X SPRI clean-ups without adapter dimer carry-over concerns.

PCR Cycles

The number of PCR cycles can be optimized to fit your particular yield requirements. In general, the cycle recommendations for total RNA can be followed when using the ribodepletion workflow. The hybridization capture workflow typically requires 200 ng of each library; thus, 4 PCR cycles should be added to the mRNA recommendations (even though total RNA is used as the input) (see **Table 5**).

Excess PCR cycles, or overamplification of the libraries, can result in the formation of a heteroduplex. This can be evidenced as a second peak >1000 bp. Heteroduplex formation can result in inaccurate quantification of library molecules and negatively impact final data quality. If a heteroduplex is observed, reduce the number of PCR cycles.

Step	Ribodepletion	Hybridization Capture
Fragmentation Time	10 min @ 94 C (Intact RNA) 5 min @ 65 C (FFPE)	15 min @ 94 C (Intact RNA) 2 min @ 94 C (FFPE)
SPRI Clean-ups	1.2X Post-RT All others the same	1.8X for ALL
PCR Cycles	Follow Total RNA recommendation	Add 4 cycles to mRNA recommendation

Table 5: Swift RNA Library Kit Protocol Adjustments for Ribodepletion or Hybridization Capture

CONCLUSION

The Swift RNA Library Kit is a fast, consistent, flexible solution for converting FFPE RNA into high-quality NGS libraries. We presented workflow options and summarized the resulting tradeoffs to data and costs, as well as highlighted protocol modifications required to implement each workflow. A compelling price per sample savings advantage is presented through the hybridization capture workflow relative to ribodepletion when exonic transcript information is the desired data output. Due to the high variability in quality and quantity of RNA derived from FFPE, this application note presents data to support the robustness of the Swift RNA Library Kit with challenging and low input RNA samples.

REFERENCES

1. Marczyk et al. 2019. The impact of RNA extraction method on accurate RNA sequencing from formalin-fixed paraffin-embedded tissues. *BMC Cancer*. 19:1189. <https://link.springer.com/article/10.1186/s12885-019-6363-0>
2. Jiang et al. 2011. Synthetic spike-in standards for RNA-Seq experiments. *Genome Research*. 21:1543-1551. <https://genome.cshlp.org/content/21/9/1543>
3. Sinicropi et al. 2012. Whole transcriptome RNA-Seq analysis of breast cancer recurrence risk using formalin-fixed paraffin-embedded tumor tissue. *PLoS ONE*. 7(7):e40092. <https://doi.org/10.1371/journal.pone.0040092>
4. RiboCop rRNA Depletion Kit V1.3. Lexogen GmbH. <https://www.lexogen.com/store/ribocop-rna-depletion-kit/>
5. Swift RNA Library Kit. Swift Biosciences. <https://swiftbiosci.com/swift-rna-library-kit/>
6. Dobin et al. 2013. STAR: ultrafast universal RNA-Seq aligner, *Bioinformatics* 29(1): 15–21. <https://doi.org/10.1093/bioinformatics/bts635>.
7. Picard tools. <http://broadinstitute.github.io/picard/>
8. DeLuca et al. 2012. RNA-SeQC: RNA-Seq metrics for quality control and process optimization. *Bioinformatics*. 28(11): 1530–1532. <https://doi.org/10.1093/bioinformatics/bts196>
9. Fgbio tools. <https://github.com/fulcrumgenomics/fgbio>.
10. Swift Exome Hyb Panel. Swift Biosciences. <https://swiftbiosci.com/swift-hybridization-capture-kits/>



Swift Biosciences, Inc.

674 S. Wagner Road, Suite 100 • Ann Arbor, MI 48103 • 734.330.2568 • www.swiftbiosci.com