

LentiStarter 2.0 kit

Components:

Reagent	Amount
pPACKH1 HIV lentiviral packaging plasmids	100 μ l
PEG-It virus precipitation reagent	20 ml
PureFection Transfection Reagent	120 μ l
TransDux virus transduction reagent	50 μ l

There will be enough reagents provided in the Kit to produce high-quality pseudoviral particles from 2 x 150mm or 5 x 100mm plates of HEK293T/FT producer cells (cells not included).

- What you will need:**
1. Your Lentivector construct (3rd generation preferable)
 2. HEK 293T or FT cells and suitable culture media
 3. Tabletop low speed centrifuge (ex. Beckman GS-6R)

Protocol

Pseudovirus Production (For 100mm plates)

Volumes for 150mm plates are denoted in **red**

Day 1

1. Plate 3×10^6 293T/FT cells (**7-8 $\times 10^6$**) in a fresh 100mm (**150mm**) plate in 10 ml (**20 ml**) of antibiotic-free DMEM medium (DMEM+FBS+Glu).

Day 2

1. The cells should be 50 to 70% confluent at day of transfection

2. Add 0.8 ml (**1.6 ml**) of serum-free DMEM media into an Eppendorf tube. Add 2 μg (**4 μg**) of transfer plasmid and 20 μl (**45 μl**) pPACKH1-plasmid mix to the tube and mix by pipetting
3. Add 24 μl (**55 μl**) of PureFection reagent to the tube, and vortex for 10 sec.
4. Incubate mixture at room temperature for a minimum of 15 min.
5. Add mixture drop-wise to the dish, and swirl to disperse evenly throughout the plates.
6. Incubate plates in 37^oC tissue culture incubator overnight.

Day 3

1. Replace transfection media with fresh complete growth media w/antibiotics

Day 4 and 5

1. Collect the medium (which now contains pseudoviral particles) into 50-ml sterile, capped conical centrifuge tubes.
2. Centrifuge at 3000 x g for 15 minutes at room temperature to pellet cell debris.
3. Transfer the viral supernatant into new fresh tubes.
4. Add **PEG-it** at a final volume of 1:5. Example: 5 ml of PEG-it should be added to 20 ml of viral supernatant, invert 10 times to mix well. Keep everything cold from this point onwards. Store virus supernatant containing PEG-it at 4^oC overnight to 3 days.

Day 6

1. Harvest PEG-it precipitated virus by centrifuging at 4^oC at 1500 x g for 30 min. Aspirate off the supernatant and resuspend the milky-white pellet in a small volume (1/100 to 1/1000 of original volume) using cold sterile PBS or cold DMEM.
2. Freeze virus aliquots at -80^oC.

Transduction of Target Cells

TransDux protocol

Day 1

1. Plate 50,000 cells per well in a 24 well plate in culture medium.

Day 2

2. Cells should be between 50 to 70% confluent.
3. Aspirate medium from cells.
4. **Add TransDux to complete medium to a final concentration of 1X (Example; add 25µl of 200x TransDux to 5 ml of medium – transfer 0.5ml of this mixture per 24 well)**
5. Add virus to each well at different MOIs or different volumes.

Day 5

6. 72 hours post transduction, the viral genome will be integrated into the host cell genome.
7. Look at the cells for reporter expression if the viral construct has a reporter like GFP and/or begin appropriate antibiotic selection to establish stable cell line.

OPTIONAL – Virus Titering

1. Aspirate off medium. Wash each well with PBS (at this point the plate can be frozen at -80°C).
2. Add 100µl of Lysis Buffer (SBI's UltraRapid Global Titering Kit) to each well.
3. Titer virus according to protocol given in the UltraRapid Global Titering Kit (SBI Cat# LV961A-1).