

chemiSOLO

Personal Western Blot Imager





Highly Sensitive Imaging

Femtogram detection of chemiluminescent Westerns



Unique Software

Capture and download images instantly from any device.




Flexible Imaging Options

Chemiluminescent blots, colorimetric SDS PAGE gels, colorimetric blots, color marker, densitometry




Chemiluminescent Blot with Color Marker




Coomassie Stained Protein Gel



Silver Stained Protein Gel



Colorimetric Stained Blot




Densitometry



Plant Bioluminescence

Quantitative Data Acquisition

chemiSOLO's Extended Dynamic Range linearly detects wider ranges without saturation.




Scan for more information on chemiSOLO's quantitative imaging



Advanced Technology



Space-efficient Design



chemiSOLO	Specifications
Part Number	CS1000
Camera	6.29MPx back-illuminated CMOS, 16-bit, 65,536 grayscale
Peltier Cooling	-40C regulated cooling from ambient
Applications	Chemiluminescent Imaging with Visible Marker, Densitometry, Visible Gel Imaging in Color and Grayscale
Light Source	LEDs
Field of View	10 x 15 cm (4.0 x 5.9 in)
Footprint (W x H x D)	29.2 x 22.2 x 43.2 cm (11.5 x 8.75 x 17.0 in)
Weight	9.1kg (20lbs)
Certifications	CE RED, cTUVus, CB scheme, RoHS

Scan for more information on the chemiSOLO





6747 Sierra Court, Suites A-B • Dublin, CA 94568 USA Phone: 925.307.7127 • Fax: 925.905.1816 info@azurebiosystems.com • www.azurebiosystems.com

Copyright © 2022 Azure Biosystems. All rights reserved. The Azure Biosystems logo and Azure Biosystems® are trademarks of the Company. All other trademarks, service marks and trade names appearing in this brochure are the property of their respective owners.